



# MEDISEPT

Perfect disinfection



# Perfect disinfection

surfaces • tools • hands

# Index

<b>Hygiene plan.....</b>	<b>3</b>
--------------------------	----------



<b>Surface disinfection.....</b>	<b>4-13</b>
----------------------------------	-------------

- Velox® Wipes NA..... 5
- Velox® Duo Wipes Tea tonic..... 6
- MEDISEPT Dry Wipes..... 7
- Velox® Spray..... 8
- Velox® Top AF..... 9
- Quatrodes® Extra..... 10
- Quatrodes® Forte..... 11
- Velox® Foam Extra..... 12
- Quatrodes® Unit NF ..... 13



<b>Instrument disinfection.....</b>	<b>15-20</b>
-------------------------------------	--------------

- Viruton® Pre..... 16
- Viruton® Pulver..... 17
- Viruton® Extra..... 18
- Viruton® Forte..... 19
- Viruton® Bohr..... 20



<b>Hand disinfection and hand care.....</b>	<b>22-30</b>
---	--------------

- MEDISEPT Wipes for hand and surface disinfection..... 23
- Velodes® Soft..... 24
- Velodes® Silk..... 25
- Velodes® Gel..... 26
- Velodes® Scrub..... 27
- Velodes® Soap..... 27
- Velodes® Soap Foam..... 28
- Dispensers and accessories..... 28/29
- Velodes® Cream..... 30
- Velodes® Cream Aloe Vera..... 30



INTENDED USE	APPLICATION	PRODUCT	DOSAGE	WHEN?	WHO?
Surfaces  Medical equipment	Cleaning and disinfection of surfaces and medical devices made of glass, metal and rubber.	<b>Quatroles® Forte</b> <b>Quatroles® Extra</b>	0,5% in 15 min. (2,5 ml per 997,5 ml of water)	Clean the surfaces after each working day or in case of need.	
	Cleaning and disinfection of surfaces sensitive to alcohol.	<b>Velox® Foam Extra</b> <b>Velox® Wipes NA</b>	Spray/wipe the surface and wait 5 min.	In case of need.	
	Cleaning and disinfection of small surfaces or surfaces difficult to access e.g. work surfaces, armchairs, dental units.	<b>Velox® Spray</b> <b>Velox® Top AF</b> <b>Velox® Wipes</b> <b>MEDISEPT wipes for hand and surface disinfection</b>	Spray/wipe the surface and wait 30 sec.	Before and after each procedure, in case of contamination or in case of need.	
	Cleaning and disinfection of suction devices.	<b>Quatroles® Unit NF</b>	2% in 15 min. (20 ml per 980 ml of water) Rinse with water.	After every patient.	
Tools  Rotary instruments	Cleaning and disinfection of dental instruments, surgical instruments and surgical devices.	<b>Viruton® Pulver</b>	1% in 30 min. (10 g per 1 l of water)	Directly after use.	
	Cleaning and disinfection of impression trays, dental braces, dental impressions.	<b>Viruton® Extra</b> <b>Viruton® Forte</b>	0,5% in 15 min. (2,5 ml per 997,5 ml of water)	Before processing.	
	Cleaning and disinfection of rotary tools, such as drills, cutters, grinding wheels.	<b>Viruton® Bohr</b>	Immerse for 15 minutes in ready-for-use product.	Directly after use. Afterwards rinse, dry and sterilise.	
	Cleaning and disinfection of dental handpieces.	<b>Velox® Spray</b> <b>Velox® Top AF</b>	Spray/wipe the surface and wait 30 sec.	After use.	
Hand disinfection  Hand washing	Hygienic and surgical hand disinfection.	<b>Velodes® Soft</b> <b>Velodes® Gel</b> <b>Velodes® Silk</b>	Hygienic disinfection in 30 sec. (2 x 3 ml liquid) Surgical disinfection in 90 sec. (2 x 3 ml liquid)	Before and after each procedure, in case of contamination or in case of need.	
	Hygienic hand wash.	<b>Velodes® Scrub</b> <b>Velodes® Soap</b>  <b>Velodes® Soap Foam</b>	Apply on wet hands so that a gentle foam is made. Rinse with water.  Apply small amount of a foam on wet hands, wash and rinse.	In case of need.	
	Hand care.	<b>Velodes® Cream</b> <b>Velodes® Cream Aloe Vera</b>	Gently massage it into the skin of hands and the nails.	Each time you wash hands or in case of need.	

# Surface disinfection



# Surface disinfection



### Logistic details

	Box	Refill
Pcs / Box	12	12
No of boxes / pallet	80	112
Pcs / europallet	960	1344

## Velox® Wipes NA

CE 2274

### Alcohol-free wipes for disinfection of small surfaces and medical equipment

Velox® Wipes NA are recommended for disinfecting surfaces sensitive to alcohol made of plastic, acrylic glass, stainless steel, metal, aluminium, rubber, and porcelain. The wipes can be used in the medical area, gastronomy, nursing homes, sanatoriums, fitness and spa centres, hairdressing and beauty salons, kindergartens, and schools. Suitable for disinfection of medical appliances, rehabilitation equipment, treatment chairs, incubators, the heads of ultrasound machines, as well as lamps in intensive care units, and operating theatres. Free from alcohol, aldehydes and phosphates, the wipes do not discolour the disinfected surfaces. The product has a pleasant scent. It is suitable for medical devices surfaces having contact with food. The product has a positive opinion of medical equipment manufacturer: Famed in terms of material tolerance for ABS and upholstery materials.

**Directions for use:** Place the refill in the box. Take a tissue out of the box and wipe a surface to be cleaned and disinfected. Close the lid. After opening the package, use the tissues within 30 days. Do not use if a tissue is dry. The surface does not need to be rinsed with water. Before applying on particularly sensitive surfaces, test the wipes on a non-visible area. For professional use only.

Wipes	Active in 1 min. against B, F, V	Medical and non-medical surfaces	0% Alcohol-free product
-------	----------------------------------	----------------------------------	-------------------------

Effectiveness	EN 14885	Conditions	
		clean	dirty
Bacteria (including MRSA)	EN 16615	60 sec.	3 min.
	EN 13727	60 sec.	60 sec.
Fungi (C. albicans)	EN 13624	60 sec.	60 sec.
	EN 16615	60 sec.	60 sec.
Tuberculosis (M. terrae)	EN 14348	5 min.	10 min.
Enveloped viruses (HIV, HBV, HCV, Influenza, Ebola, Vaccinia, BVDV, Herpes simplex, SARS-CoV-2)	EN 14476	60 sec.	60 sec.

Active ingredients: amine, QAC  
Packaging units: 100 pieces box, 100 pieces refill

Non-woven wipes 13x20cm, 23g/m<sup>2</sup>  
Permission number for biocidal products in Poland: 8023/20

# Surface disinfection



### Logistic details

	Box	Refill
Pcs / Box	15	12
No of boxes / pallet	63	112
Pcs / europallet	945	1344

## Velox® Duo Wipes Tea tonic

CE 2274

### Alcohol wipes for cleaning and rapid disinfection of small surfaces and medical equipment

Velox® Duo Wipes are destined for rapid cleaning and disinfection of small surfaces in the medical area, nursing homes, sanatoriums, fitness and SPA centres, hairdressing, massage and beauty salons, kindergartens, and schools. The product is suitable for surfaces having contact with food. Recommended for small surfaces of medical equipment: operating beds and chairs, beds, medical and surgical equipment, and disinfection of dental ends before the sterilisation process, prosthetic spoons, and silicon moulds. Disinfect the surfaces killing in 30 seconds bacteria, fungi, mycobacteria, enveloped viruses (such as Vaccinia, BVDV, SARS-CoV-2, HIV, HBV, HCV), Rota and Noro viruses.

**Directions for use:** Take a wipe out of the package and wipe a surface to be cleaned and disinfected. Close the package tightly. After opening the package, use the wipes within 21 days. Do not use if the wipe is dry. In case of disinfecting surfaces having contact with food, rinse them with drinking water. Other surfaces do not require rinsing with water. Before applying on particularly sensitive surfaces, test the wipes on a non-visible area. Do not use on surfaces sensitive to alcohol. For professional use only.

	50 pcs packaging		Active against Tbc in 5 min.		Medical and non-medical surfaces		Tea tonic scent
--	------------------	--	------------------------------	--	----------------------------------	--	-----------------

Effectiveness	EN 14885	Conditions	
		clean	dirty
Bacteria (including MRSA)	EN 13727	30 sec.	30 sec.
	EN 16615	30 sec.	30 sec.
Fungi (C. albicans)	EN 13624	30 sec.	30 sec.
	EN 16615	30 sec.	30 sec.
Mycobacteria tuberculosis (M. terrae)	EN 14348	30 sec.	30 sec.
Enveloped viruses (HIV, HBV, HCV, Vaccinia, Herpes simplex, Ebola, SARS-CoV-2)	EN 14476	30 sec.	60 sec.
BVDV virus	EN 14476	30 sec.	30 sec.
Rota virus	EN 14476	30 sec.	60 sec.
Noro virus	EN 14476	60 sec.	30 sec.

Active ingredients: ethanol, propan-2-ol  
Packaging units: 50 pcs flow pack with self-adhesive closing

Non-woven wipes: 19x15cm, 50 g/m<sup>2</sup>  
Permission number for biocidal products in Poland: 7905/19

# Surface disinfection



## MEDISEPT Dry Wipes MEDISEPT Dry Wipes Box



Dry wipes for soaking with disinfecting/cleaning agent

MEDISEPT Dry Wipes are designed for soaking with a ready-to-use product or utility solution of concentrate with water. Formulated for cleaning and disinfecting various surfaces of non-invasive medical devices. The wipes may be used in healthcare facilities, surgeries, nursing homes, and other areas, according to the intended use of the soaking agent.

**Directions for use:** Prepare the soaking liquid according to the recommendations on the product's label. Wipe the surface with a damp wipe and leave it to dry for the time recommended by the soaking-liquid's manufacturer. Clean the bucket and lid before re-use with the soaking solution. Using the liquid-soaked wipes does not require protective gloves.

 <p>Dry wipes 100 pcs</p>	 <p>For small and large surfaces</p>	 <p>Durable and convenient bucket</p>
--	---	--

Products recommended for use with MEDISEPT Dry Wipes:



The compatibility of products with the wipes has been confirmed by microbiological-efficiency test using the four-field test according to the EN 16615 standard.

Composition of nonwovens: polypropylene (50%), cellulose (50%)  
 Packaging units:  
 - bucket with wipes for soaking insert, 100 pcs., EAN 5907626638035  
 - supplementary insert for the bucket, 100 pcs., EAN 5907626638059

Wipe grammage: 45 g/m<sup>2</sup>  
 Wipe size: 16x30cm  
 Bucket material: HDPE with high-resistance to chemicals

# Surface disinfection



### Logistic details

	1 L	5 L
Pcs / Box	12	1
No of boxes / pallet	32	96
Pcs / europallet	384	96

## Velox® Spray

CE 2274

### Spray for surfaces cleaning and disinfection

Velox® Spray is an alcoholic liquid for rapid disinfection of small surfaces in the medical area, nursing homes, sanatoriums, fitness and SPA centres, hairdressers, massage and beauty salons, kindergartens, and schools. The product is suitable for surfaces having and not-having contact with food. It can be used for disinfecting surfaces such as table tops, tables, handles, for hairdressing and cosmetics equipment and accessories. The product is recommended also for disinfection of non-invasive medical devices such as treatment chairs and medical equipment. Velox® Spray is recommended for the disinfection of dental handpieces (before sterilisation), prosthetic trays and silicone dental impressions. Velox® Spray has bactericidal, yeasticidal, cidal against mycobacteria, and limited virucidal activity (against Rota virus, Noro virus) and against enveloped viruses (including Vaccinia, BVDV, HIV, HBV, HCV, SARS-CoV-2).

**Directions for use:** Spray the surface from the distance of about 30 cm, ensuring it is fully covered. Leave it for 30 seconds. In case of disinfecting surfaces having contact with food, rinse them with drinking water. Other surfaces do not require rinsing with water. Do not use on surfaces sensitive to alcohol (e.g. acrylic glass). For professional use only.

Spray	Hairdressing and cosmetics tools	Neutral scent	Tea tonic scent
-------	----------------------------------	---------------	-----------------

Effectiveness	EN 14885	Conditions	
		clean	dirty
Bacteria (including MRSA)	EN 13727	30 sec.	30 sec.
	EN 16615	30 sec.	30 sec.
Fungi (C. albicans)	EN 13624	15 sec.	15 sec.
	EN 16615	30 sec.	30 sec.
Mycobacteria tuberculosis (M. terrae)	EN 14348	30 sec.	30 sec.
Enveloped viruses (Vaccinia, HIV, HBV, HCV, Herpes simplex, Ebola, SARS-CoV-2)	EN 14476	30 sec.	60 sec.
BVDV	EN 14476	30 sec.	30 sec.
Rota virus	EN 14476	30 sec.	60 sec.
Noro virus	EN 14476	30 sec.	30 sec.

Active ingredients: ethanol, propan-2-ol  
Packaging units: 1L spray bottle, 5L canister

Permission number for biocidal products in Poland: 7675/19

# Surface disinfection



### Logistic details

	1 L	5 L
Pcs / Box	12	1
No of boxes / pallet	32	96
Pcs / europallet	384	96

## Velox® Top AF

CE 2274

### Spray for surfaces cleaning and disinfection

Velox® Top AF is an alcohol based formula designed for rapid disinfection of small surfaces such as operating chairs and tables, medical equipment, and dental tips before sterilization. It is designed for usage in gastronomy, nursing homes, sanatoriums, gyms, fitness and SPA centres, hairdressing, massage and beauty salons, tattoos parlours, kindergartens, and schools. The product is suitable for surfaces having contact with food. Ideal for cleaning small and hard-to-reach surfaces such as tops, tables, handles, for hairdressing and cosmetic equipment and accessories. The product has a bactericidal, yeasticidal, cidal against mycobacteria, and limited virucidal activity. Active against Rota and Adeno viruses. The product has a positive opinion of medical equipment manufacturer: Famed in terms of material tolerance for ABS and upholstery materials.

**Directions for use:** Spray the surface from the distance of 30 cm covering the whole disinfected area. Spread the liquid using disposable towels and allow to dry. The surface does not need rinsing with water, dries out completely and does not leave smudges. Do not use on surfaces sensitive to alcohol (i.e. acrylic glass). For professional use only.

Spray	Active against Tbc, Rota, Adeno	Neutral scent	Grapefruit scent
-------	---------------------------------	---------------	------------------

Effectiveness	EN 14885	Conditions	
		clean	dirty
Bacteria (including MRSA)	EN 1276	30 sec.	60 sec.
	EN 13697	30 sec.	60 sec.
	EN 13727	30 sec.	
Fungi (C. albicans)	EN 1650	30 sec.	60 sec.
	EN 13624	30 sec.	
Tuberculosis (M. terrae)	EN 14348	30 sec.	60 sec.
Enveloped viruses (HIV, HBV, HCV, Vaccinia, Influenza, Ebola, BVDV, SARS-CoV-2)	DVV/RKI	30 sec.	30 sec.
Rota virus	EN 14476	30 sec.	
Adeno virus	EN 14476	60 sec.	

Active ingredients: propan-2-ol, ethanol, amine, QAC  
Packaging units: 1L spray bottle, 5L canister

Permission number for biocidal products in Poland: 6140/15

# Surface disinfection



### Logistic details

	1 L	5 L
Pcs / Box	12	1
No of boxes / pallet	32	96
Pcs / europallet	384	96

## Quatrodes® Extra

CE 2274

### Concentrate for surfaces cleaning and disinfection

Quatrodes® Extra is a liquid concentrate for simultaneous cleaning and disinfection of surfaces in the medical area (units, operating chairs and beds, operating tables, medical equipment, floors, and walls), in gastronomy, nursing homes, sanatoriums, gyms, fitness and SPA centres, hairdressing, massage and beauty salons, solariums (for disinfecting sunbeds), tattoos parlours, kindergartens, and schools. The product is suitable for surfaces having contact with food. The substance can be used on all floor and above floor surfaces (doors, tabletops, cabinets, tables) and for cleaning and disinfecting hairdressing accessories.

**Directions for use:** Prepare a working solution by diluting the concentrate with water in a desired proportion. The choice of concentration depends on the desired spectrum and duration of operation. In order to obtain 0,5% concentration, dilute 5 ml of the concentrate in 995 ml of water. Disinfected surfaces do not require rinsing with water. When disinfecting surfaces are to be made contact with food, rinse them with drinking water. Do not combine with other detergents nor disinfectants. For professional use only.

Concentrate	Active against B, F, V, Tbc, Adeno, Polio	Medical and non-medical surfaces	Concentration from 0,5% in 15 min.
-------------	---	----------------------------------	------------------------------------

Effectiveness	EN 14885	Time	
		15 min.	30 min.
Bacteria (including MRSA)	EN 13727	0,25%*, 0,5%	
	EN 14561	1%	0,5%
Fungi (C. albicans)	EN 13624	0,25%*, 0,5%	
	EN 14562	1%	0,5%
Mycobacteria (M.terrae, M. avium)	EN 14348	0,5%*, 1%	0,5%
	EN 14563	1%*, 2%	0,5%
Enveloped viruses (including HBV, HIV, HCV, Vaccinia, BVDV, Herpes simplex, influenza A, B, C, Ebola, SARS-CoV-2)	EN 14476	0,5%	
Adeno virus	EN 14476	0,5%*	0,5%
Polio virus	EN 14476	1%	0,5%*

\*clean conditions

Active ingredients: amine, QAC  
Packaging units: 1L bottle with doser, 5L canister

Dosage: for 0,5% -> 5ml in 995ml of water  
Permission number for biocidal products in Poland: 6046/15

# Surface disinfection



### Logistic details

	1 L	5 L
Pcs / Box	12	1
No of boxes / pallet	32	96
Pcs / europallet	384	96

## Quatrodes® Forte

CE 2274

### Concentrate for surfaces cleaning and disinfection

Quatrodes® Forte is a liquid concentrate for simultaneous cleaning and disinfection of surfaces in the medical area (units, operating chairs and beds, operating tables, medical equipment, floors, and walls), in gastronomy, nursing homes, sanatoriums, gyms, fitness and SPA centres, hairdressing, massage and beauty salons, solariums (for disinfecting sunbeds), tattoos parlours, kindergartens, and schools. The product is suitable for surfaces having contact with food. The product has a bactericidal, mycobactericidal, yeast- and virucidal effect. The substance can be used on all floor and above floor surfaces (doors, tabletops, cabinets, tables), as well as for cleaning and disinfecting hairdressing accessories.

**Directions for use:** Prepare a working solution by diluting the concentrate with water in a desired proportion. The choice of concentration depends on the desired spectrum and duration of operation. In order to obtain a concentration of 1% dilute 10 ml of the concentrate in 990 ml of water. Disinfected surfaces do not require rinsing with water. Do not combine with other detergents nor disinfectants. For professional use only.

	Concentrate		Active against B, F, V, Tbc, Adeno, Polio		Medical and non-medical surfaces		Concentration from 0,5% in 15 min.
--	-------------	--	---	--	----------------------------------	--	------------------------------------

Effectiveness	EN 14885	Time	
		15 min.	30 min.
Bacteria (including MRSA)	EN 13727	0,5%*, 1%	
	EN 14561	1%	0,5%
	EN 13697	1%*	
Fungi (C. albicans)	EN 13624	0,5%	
	EN 14562	1%	0,5%
	EN 13697	1%*	
Mycobacteria (M. terrae, M. avium)	EN 14348, EN 14563	4%	1%
Enveloped viruses (Vaccinia, BVDV, HBV, HIV, HCV, Herpes simplex, Influenza A, Ebola, SARS-CoV-2)	EN 14476	0,5%	
Adeno virus	EN 14476	1%	
Polio virus	EN 14476	4%	

\*clean conditions

Active ingredients: amine, QAC  
Packaging units: 1L bottle with doser, 5L canister

Dosage: for 0,5% -> 5ml in 995ml of water  
Permission number for biocidal products in Poland: 6066/15

# Surface disinfection



### Logistic details

	1 L	5 L
Pcs / Box	12	1
No of boxes / pallet	32	96
Pcs / europallet	384	96

## Velox® Foam Extra

CE 2274

### Foam for cleaning and rapid disinfection of sensitive surfaces

Velox® Foam Extra is recommended for cleaning and disinfection of a wide range of surfaces in medical area, nursing homes, fitness and spa centres, massage and beauty salons. It may be used in hospitals at intensive care units and surgery theatres for disinfection of medical devices, rehabilitation equipment, treatment chairs, incubators, ultrasound heads, and lamps. Free from aldehydes and phosphates, does not discolour disinfected surfaces. The product has a pleasant, fresh scent. It is suitable for surfaces having contact with food. The product has a positive opinion of medical equipment manufacturer: Famed in terms of material tolerance for ABS and upholstery materials.

**Directions for use:** The product has a double spray and foam head. Application in foam is recommended. Spray the surface from a distance of 20-30 cm, making sure it is completely covered. Remove the excess foam with a disposable paper towel. The surface does not require rinsing with water. Before disinfecting particularly sensitive surfaces (e.g. natural leather) it is recommended to perform the test on a non - visible part of the surface. For professional use only.

Foam	Active against Tbc in 5 min.	Medical and non-medical surfaces	0% Alcohol-free product
------	------------------------------	----------------------------------	-------------------------

Effectiveness	EN 14885	Conditions	
		clean	dirty
Bacteria (including MRSA)	EN 16615	60 sec.	3 min.
	EN 13727	60 sec.	60 sec.
Fungi ( <i>C. albicans</i> )	EN 13624	60 sec.	60 sec.
	EN 16615	60 sec.	60 sec.
Tuberculosis ( <i>M. terrae</i> )	EN 14348	5 min.	10 min.
Enveloped viruses (HIV, HBV, HCV, Influenza, Ebola, Vaccinia, BVDV, Herpes simplex, SARS-CoV-2)	EN 14476	60 sec.	60 sec.

Active ingredients: amine, QAC  
Packaging units: 1L spray bottle, 5L canister

Permission number for biocidal products in Poland: 8192/20

# Surface disinfection



### Logistic details

	1 L	5 L
Pcs / Box	12	1
No of boxes / pallet	32	96
Pcs / europallet	384	96

## Quatrodes® Unit NF

CE 2274

### Concentrate for cleaning and disinfection of suction devices

Quatrodes® Unit NF is designed for simultaneous cleaning and disinfection of suction systems, amalgam separators, and dental unit bowls. The product does not foam. Used daily, it prevents biofilm from forming. It can also be used to disinfect closed circuit in hydromassage tubs that are medical devices. The solution can remain in circulation overnight. Before the first patient, it is recommended to rinse the closed system with water.

**Directions for use:** The concentrate should be diluted with water. In order to obtain 2% concentration, dilute 20 ml of the concentrate in 980 ml of water. **In suction system:** Introduce the solution to the suction system, depending on the degree of soiling, once or twice a day. Suck the whole product, turn the device off and leave it for at least 15 minutes. After disinfection, rinse the suction system with water. Rinsing can be carried out after the recommended disinfection time or the next day, immediately after turning the unit on. Do not combine with other detergents nor disinfectants. **In hydromassage tub:** Before the first use, rinse the device with water to remove any possible residues of other product. Fill the tub with clean, cold water to the level above the nozzles. Add 200 ml of the product for every 10 L of water. Turn the hydromassage system on for 15 minutes. After use, empty the tub and rinse with water. Do not combine with other agents. Do not use the product during the patient's healing bath.

Concentrate	Active against B, F, V, Tbc	Dental unit	Concentration from 2% in 15 min.
-------------	-----------------------------	-------------	----------------------------------

Effectiveness	EN 14885	Time		
		15 min.	30 min.	60 min.
Bacteria (including MRSA)	EN 13727	1%		
Fungi (C. albicans)	EN 13624	1%*		
	EN 13624	2%	1%	
Tuberculosis (M. terrae)	EN 14348			4%
Enveloped viruses (HBV, HIV, HCV, Ebola, Vaccinia, Influenza, Herpes simplex, SARS-CoV-2)	EN 14476	2%		

\* clean conditions

Active ingredients: amine, QAC

Dosage: for 1% -> 10ml in 990ml of water

Packaging units: 1L bottle with doser, 5L canister

## The main features of the packaging

**Integrated dispenser for the easy dosing of a concentrate**

**High-quality white plastic**

**The background colour of the label corresponding to the application**

**Sections for product-form Pictograms, target areas and additional information**

**QR code on the back which, after scanning, redirects you to the product's website.**

**Product name**

**Brief description of the product**

**Dosage instructions**

**Product-category code**

## Colour-coding of products

Target areas:

Surfaces	Instruments	Hands
----------	-------------	-------

Intended use category:

Disinfection	Care / Washing
--------------	----------------

# Instrument disinfection



# Instrument disinfection



### Logistic details

	1 L
Pcs / Box	12
No of boxes / pallet	32
Pcs / europallet	384

## Viruton® Pre

CE 2274

### Pre-disinfecting and cleaning foam for surgical instruments

Viruton® Pre is a ready-to-use foam for pre-disinfection and cleaning of instruments before the actual disinfection process. The product contains corrosion inhibitors and has a high material tolerance making it ideal for instruments made of stainless steel, galvanized steel, aluminium, plastic, and rubber. The product has very good cleaning properties also in case of dried organic residues. Viruton® Pre prevents residues from drying during collection and transportation of instruments to a place of proper disinfection. It has a pleasant scent, and both bacteriostatic and bactericidal properties. Viruton® Pre has a broad biocidal spectrum and exceptional cleaning properties thanks to the innovative formula based on the synergy of surfactants and nano-micelle technology.

**Directions for use:** Place the tools immediately after use in the container and spray thoroughly with the undiluted product, making sure that all elements are covered with foam. After a certain period of time, remove the tools from the tray and rinse with water. Pre-disinfection process may take up to 48 hours. In case of visible contamination, clean the instruments with a brush and cover again with the foam. Once rinsed with water, instruments shall be disinfected either manually or automatically. Do not combine with other detergents nor disinfectants. Use protective gloves. Do not use on instruments made of copper, brass nor zinc. For professional use only.

Foam	Active against B, F, V, Tbc, Adeno, Polio	Instruments	Active in 15 min.
------	---	-------------	-------------------

Effectiveness	EN 14885	Dirty conditions
Bacteria (including MRSA)	EN 13727	15 min.
	EN 14561	15 min.
Fungi (C. albicans)	EN 13624	15 min.
	EN 14562	15 min.
Mycobacteria (M. terrae, M. avium)	EN 14348	15 min.
	EN 14563	15 min.
Enveloped viruses (Vaccinia, BVDV, HBV, HIV, HCV, Ebola, SARS-CoV-2)	EN 14476	15 min.
Adeno virus	EN 14476	15 min.
Polio virus	EN 14476	15 min.

Active ingredients: amine, QAC, corrosion inhibitors

Packaging units: 1 L bottle with DUO spray

# Instrument disinfection



### Logistic details

	1 kg	5 kg
Pcs / Box	1	1
No of boxes / pallet	240	80
Pcs / europallet	240	80

## Viruton® Pulver

CE 2274

### Powder for cleaning and disinfecting tools, endoscopes, incubators and other surfaces of medical devices

Viruton® Pulver is a modern generation preparation for cleaning and high-level disinfection of invasive and non-invasive medical devices. It is recommended for surgical and dental instruments, all kinds of endoscopes, surgical probes, ventilator tubes, and other anesthetic equipment. It can also be used for the disinfection of surfaces of medical equipment, such as incubators. It is recommended for manual disinfection, in ultrasonic cleaners, automatic washers, and endoscope cleaning baths. The product does not require an activator. Untampered solution is stable up to 30 hours. The innovative formula contains four enzymes (protease, lipase, amylase and mannanase) responsible for removing organic contaminants, such as blood, protein, fats, and secretions. The composition of surfactants prevents the powder from dusting and guarantees excellent washing properties. The product shows anti-corrosion properties and high material tolerance - opinion of tool manufacturer, Chirmed.

**Directions for use:** Refer to the dosing table and measure the required amount of tap water. Add the correct amount of powder, using the provided measuring cup, stir the solution. Cover the container and leave it for 15 minutes to activate. Afterwards, stir the solution again. Using warm water will speed up the dissolution of the powder. Small amounts of residue on the bottom of the container ensure the stability of the solution in time. Change the solution daily. Test strips available to confirm the activity of the working solution. Detailed instructions for disinfection are available.

Non - dusting	Active against Polio, Adeno, Noro, Spores	Active washing	Concentration from 1% in 30 min.
---------------	---	----------------	----------------------------------

Effectiveness	EN 14885	Time		
		10 min.	30 min.	60 min.
Bacteria (S. aureus, P. aeruginosa, E. hirae, MRSA)	EN 13727		0,5%	
	EN 14561	1%	0,5%	
	EN 13697		0,5%	
Fungi (C. albicans)	EN 13624		0,5%	
	EN 14562	1%	0,5%	
Fungi (A. Brasiliensis)	EN 13624	1%		0,5%
	EN 14562	2%		
Tuberculosis (M. terrae)	EN 14348		0,5%*	
	EN 14563			0,5%
Mycobacteria (M. terrae, M. avium)	EN 14476	1%	1%	
Viruses (Polio, Adeno, Noro, HBV, HIV, HCV, Rota, Influenza, Herpes simplex, Ebola, SARS-CoV-2)	EN 14476	2%	1%	
Spores (C. difficile)	EN 13704	1%*, 2%	0,5%	
Spores (B. subtilis)	EN 13704	1%*	0,5%	

\*clean conditions

Active ingredients: sodium percarbonate, TAED, enzymes, organic acids, non-ionic surfactants, corrosion inhibitors

Packaging units: 1kg with measuring cup, 5kg with measuring cup (use the measuring cup according to the dosing table)

# Instrument disinfection



### Logistic details

	1 L	5 L
Pcs / Box	12	1
No of boxes / pallet	32	96
Pcs / europallet	384	96

## Viruton® Extra

CE 2274

### Concentrate for cleaning and disinfection of surgical instruments

Viruton® Extra is a highly effective concentrate for simultaneous cleaning and disinfecting surgical and rotary instruments. It contains an enzyme – protease, which provides excellent cleaning properties. Designed for manual disinfection, as well as in ultrasonic cleaners. Effectively removes dried out organic soiling. It is recommended for instruments made of stainless steel, nickel, copper, aluminium, rubber, porcelain, glass, and plastic. Thanks to an innovative formula based on the synergy of surfactants and nano-micelli technology, the product has a wide range of effectiveness and exceptional cleaning properties. The product shows anti-corrosion properties and high material tolerance - opinion of tool manufacturer, Chirmed. For hygienic reasons, a daily change of the solution is recommended.

**Directions for use:** Prepare working solution by diluting the concentrate with water to a desired concentration. The level of dilution depends on desired spectrum and duration of operation. In order to obtain 0.5% concentration, dilute 5 ml of the concentrate with 995 ml of water. Immerse instruments in the solution immediately after use. All parts and holes must be completely covered by the solution. Once disinfected, rinse the instruments with water and dry. The uncontaminated working solution is active for 14 days. For hygiene reasons, it is recommended to exchange the working solution on a daily basis. Do not combine with other detergents or disinfectants. For professional use only.

Concentrate	Active against Adeno, Polio	Contains protease enzyme	Concentration from 0,5% in 15 min.
-------------	-----------------------------	--------------------------	------------------------------------

Effectiveness	EN 14885	Time	
		15 min.	30 min.
Bacteria (including MRSA)	EN 13727	0,5%	
	EN 14561	1%	0,5%
Fungi (C. albicans)	EN 13624	0,5%	
	EN 14562	1%	0,5%
Mycobacteria (M. terrae, M. avium)	EN 14348	0,5%*, 1%	0,5%
	EN 14563	1%*, 2%	0,5%
Enveloped viruses (Vaccinia, BVDV, HBV, HIV, HCV, Herpes simplex, Influenza type A, Ebola, SARS-CoV-2)	EN 14476	0,5%	
Adeno virus	EN 14476	0,5%*	0,5%
Polio virus	EN 14476	1%	0,5%*

\*clean conditions

Active ingredients: amine, QAC, enzyme, corrosion inhibitors  
Packaging units: 1L bottle with doser, 5L canister

Dosage: for 0,5% -> 5ml in 995ml of water

# Instrument disinfection



### Logistic details

	1 L	5 L
Pcs / Box	12	1
No of boxes / pallet	32	96
Pcs / europallet	384	96

## Viruton® Forte

CE 2274

### Concentrate for cleaning and disinfection of surgical instruments

Viruton® Forte is a highly effective concentrate for simultaneous cleaning and disinfection of surgical and rotary instruments. It has a very good cleaning properties and wide microbiological effectiveness within a short period of time. Designed for manual instrument disinfection, as well as in ultrasonic cleaners. The product presents anti-corrosion properties and a high material tolerance. Viruton® Forte is recommended for instruments made of stainless steel, nickel, copper, aluminium, rubber, porcelain, glass, and plastics. Thanks to an innovative formula based on the synergy of surfactants and nano-micelli technology, the product has a wide range of effectiveness and exceptional cleaning properties. The product shows anti-corrosion properties and high material tolerance - opinion of tool manufacturer, Chirmed. For hygienic reasons, a daily change of the solution is recommended.

**Directions for use:** Prepare working solution by diluting the concentrate with water. The level of dilution depends on a desired spectrum and duration of operation. In order to obtain a concentration of 1% dilute 10 ml of the concentrate in 990 ml of water. Instruments shall be immersed in the solution immediately after use. All parts and holes must be completely filled in and covered by the solution. Once disinfected, rinse the instruments with water and dry. The uncontaminated working solution is active for 14 days. For hygiene reasons, it is recommended to exchange the working solution on a daily basis. Do not combine with other detergents or disinfectants. For professional use only.

Concentrate	Active against B, F, V, Tbc, Adeno, Polio	Instruments	Concentration from 0,5% in 15 min.
-------------	---	-------------	------------------------------------

Effectiveness	EN 14885	Time	
		15 min.	30 min.
Bacteria (including MRSA)	EN 13727	0,5%*, 1%	
	EN 14561	1%	0,5%
	EN 13697	1%*	
Fungi (C. albicans)	EN 13624	0,5%	
	EN 14562	1%	0,5%
	EN 13697	1%*	
Mycobacteria (M. terrae, M. avium)	EN 14348, EN 14563	4%	1%
Enveloped viruses (Vaccinia, BVDV, HBV, HIV, HCV, Herpes simplex, Influenza type A, Ebola, SARS-CoV-2)	EN 14476	0,5%	
Adeno virus	EN 14476	1%	
Polio virus	EN 14476	4%	*clean conditions

Active ingredients: amine, QAC, corrosion inhibitors  
Packaging units: 1L bottle with doser, 5L canister

Dosage: for 1% -> 10ml in 990ml of water

# Instrument disinfection



### Logistic details

	1 L	5 L
Pcs / Box	12	1
No of boxes / pallet	32	96
Pcs / europallet	384	96

## Viruton® Bohr

CE 2274

### Ready-to-use formula for cleaning and disinfection of rotary instruments

Viruton® Bohr is recommended for simultaneous cleaning, and disinfection of precise rotary instruments and other invasive medical tools. Designed for drills and milling cutters made of steel, hard metals, diamonds, abrasive and polishing tools, such as silicone rubbers, ceramic tiles, and root canal treatment instruments. Viruton® Bohr contains inhibitors of corrosion, which protect instruments made of steel and metal alloys. It has excellent cleaning properties also in case of dried organic residues. Instruments can be disinfected by immersion in a disinfection tub or in an ultrasonic cleaner.

**Directions for use:** Place the instruments in the undiluted formula immediately after use. Carefully immerse all elements. After given period of time, take the instruments out, rinse with water and dry. In case of visible contamination, clean it with a brush and immerse the tools in the formula again. After drying, sterilize the instruments. Do not combine the product with other detergents nor disinfectants. For professional use only.

Liquid	Active against B, F, V, Tbc, Adeno, Polio	Rotary instruments	Active in 15 min.
--------	---	--------------------	-------------------

Effectiveness	EN 14885	Dirty conditions
Bacteria (including MRSA)	EN 13727	15 min.
	EN 14561	15 min.
Fungi (C. albicans)	EN 13624	15 min.
	EN 14562	15 min.
Mycobacteria (M. terrae, M. avium)	EN 14348	15 min.
	EN 14563	15 min.
Enveloped viruses (Vaccinia, BVDV, HBV, HIV, HCV, Herpes simplex, Influenza type A, Ebola, SARS-CoV-2)	EN 14476	15 min.
Adeno virus	EN 14476	15 min.
Polio virus	EN 14476	15 min.

Active ingredients: amine, QAC, corrosion inhibitors

Packaging units: 1L bottle with doser, 5L canister

# How to prepare a working solution?

In order to obtain selected concentration, measure the appropriate amount of water in a tub or a bucket, then add the amount of disinfectant concentrate given in the table.

Concentration	4 %	2 %	1 %	0,5 %	0,25 %
Concentrate amount (ml)	40	20	10	5	2,5
Water amount (ml)	960	980	990	995	997,5

## Concentrate efficiency (liters of solution)

Concentration	Efficiency from 1l of concentrate	Concentration	Efficiency from 1l of concentrate
0,25 %	400 liters	2 %	50 liters
0,5 %	200 liters	3 %	33 liters
1 %	100 liters	4 %	25 liters

## How to use a doser?

**1** Open the cap over the dispenser

**2** Press the bottle. The liquid flows through the tube to the appropriate level of the dispenser

**3** Tilt the bottle as shown in the picture and pour the discharged amount of concentrate into the water



# Hand disinfection and hand care



# Hand disinfection and hand care



## Logistic details

	flow pack
Pcs / Box	15
No of boxes / pallet	63
Pcs / europallet	945

## MEDISEPT Wipes for hand and surface disinfection

High-quality wipes for disinfecting hands and all small surfaces

MEDISEPT Wipes are high-quality alcoholic wipes designed for hygienic hand disinfection and any small alcohol-resistant surfaces. They are recommended for surfaces such as tabletops, knobs, handles, seating surfaces, medical equipment, and accessories made of plastic, stainless steel, porcelain, and ABS plastic. Due to its intended use on the skin, the product contains care-substance additives, which moisturize the hands. It has bactericidal, yeasticidal, mycobactericidal, tuberculocidal, and limited virucidal effect (e.g. Herpes Simplex, HIV, Hepatitis B and C, Rota virus). Dermatologically tested.

**Directions for use:** Take a wipe out of the package and wipe hands or surface. Close the package tightly. After opening, use the wipes within 21 days. Do not use if the wipe is dry. Do not use on surfaces sensitive to alcohol.

	50 pcs packaging		Active against Tbc, Rota in 30 sec.		Disinfect and protect hands		Disinfects surfaces
--	------------------	--	-------------------------------------	--	-----------------------------	--	---------------------

Effectiveness	EN 14885	Time
Bacteria (including MRSA)	EN 13727	30 sec.
Bacteria	EN 16615	30 sec.
Fungi (C. albicans)	EN 1275, EN 13624	30 sec.
M. Tuberculosis (M. terrae)	EN 16615	30 sec.
M. Tuberculosis (M. terrae)	EN 14348	30 sec.
Enveloped viruses (HIV, HBV, HCV, Influenza, Vaccinia, BVDV, Herpes simplex, Ebola, SARS-CoV-2)	EN 14476	30 sec.
Rota Virus	EN 14476	45 sec.
Hygienic disinfection*	EN 1500	30 sec.

\*test performed on saturated formulation, 1 tissue contains approx. 3.56 ml of disinfecting liquid

Active ingredients: propan-2-ol, d-gluconic acid, compound with N, N'-bis(4-chlorophenyl)-3,12-diimino-2,4,11,13-tetraazatetradecanediamidine (2:1) (chlorhexidine gluconate)

Packaging units: 50 pcs flow pack with self-adhesive closing

Non-woven fabric: 19x15cm, 50 g/m<sup>2</sup>

Permission number for biocidal products in Poland: 6678/16

# Hand disinfection and hand care



### Logistic details

	250 ml	500 ml	1 L	1 L bag	5 L
Pcs / Box	40	24	12	6	1
No of boxes / pallet	30	40	40	80	96
Pcs / europallet	1200	960	480	480	96

## Velodes® Soft

### Liquid for hygienic and surgical hand disinfection

Velodes® Soft is a ready-to-use, alcohol based product designed for hygienic and surgical disinfection of hands. It has a wide and confirmed disinfection properties against bacteria, mycobacteria, fungi, and some viruses. It has prolonged bactericidal effect. It deactivates enveloped viruses such as HIV, HBV, HCV, SARS-CoV-2 and Rota virus. Contains glycerine, which moisturizes and cares for hands. The special formula prevents from hand stickiness. Dermatologically tested.

**Directions for use:** Hygienic hand disinfection (according to EN 1500): rub 2 x 3 ml of the product in hands for 30 seconds. Surgical hand disinfection (according to EN 12791): rub 2 x 3 ml of the product in clean and dry hands for 90 seconds. Keep the hands equally moist by adding more liquid. The product has a prolonged surgical disinfection effect up to 3 hours. For professional use only.

	Liquid		Hygienic disinfection in 30 sec.		Surgical disinfection in 90 sec.		Prolonged effect
--	--------	--	----------------------------------	--	----------------------------------	--	------------------

Effectiveness	EN 14885	Time
Bacteria (including MRSA)	EN 1500	30 sec.
	EN 12791	90 sec.
	EN 13727	30 sec.
Fungi (C. albicans)	EN 1275	30 sec.
	EN 13624	30 sec.
Tuberculosis (M.terrae)	EN 14348	30 sec.
Enveloped viruses (Vaccinia and BVDV and other like: HBV, HIV, HCV, SARS-CoV-2, Herpes, Influenza type A, Ebola)	EN 14476	30 sec.
Rota virus	EN 14476	45 sec.

Active ingredients: propan-2-ol, chlorhexidine gluconate  
 Packaging units: 250ml spray bottle, 500ml bottle, 1L bottle, 1L bag, 5L canister

Permission number for the biocidal product in Poland: 6174/15

# Hand disinfection and hand care



### Logistic details

	500 ml	1L bag	5L
Pcs / Box	24	12	1
No of boxes / pallet	40	40	96
Pcs / europallet	960	480	96

## Velodes® Silk

### Liquid for hygienic and surgical hand disinfection, as well as surface disinfection

Velodes® Silk is a ready-to-use, alcohol-based liquid for hygienic and surgical hand disinfection in the medical field, the food sector, industrial, domestic environments, and in public utilities. For surfaces having and not-having contact with food. Designed for disinfection of non-invasive medical devices. The product has a bactericidal, yeasticidal, and limited virucidal activity (against Rota and Noro viruses), it also works against tuberculosis and enveloped viruses (including Vaccinia, BVDV, HIV, HBV, HCV, Influenza, and SARS-CoV-2). Liquid form.

**Directions for use:** Hygienic hand disinfection (according to EN 1500): rub 2 x 3 ml of the product to for 30 seconds (2 x 15 sec.). Surgical hand disinfection (according to EN 12791): rub 2 x 3 ml of the product for 3 minutes (2 x 90 sec.). Keep hands equally moist by adding more liquid. Surface disinfection: Spray the surface from the distance of approximately 30 cm and leave it to dry. Do not apply on alcohol-sensitive surfaces (e.g. acrylic glass). For professional use only.

	Liquid		Active against B, F, V, Tbc, Rota, Noro		Hygienic hand disinfection in 30 sec.		Surgical hand disinfection in 3 min.
--	--------	--	---	--	---------------------------------------	--	--------------------------------------

Effectiveness	EN 14885	Conditions	
		clean	dirty
Bacteria (including MRSA)	EN 13727	30 sec.	30 sec.
Fungi (C. albicans)	EN 13624	15 sec.	15 sec.
Tuberculosis (M. terrae)	EN 14348	30 sec.	30 sec.
Enveloped viruses (Vaccinia, HIV, HBV, HCV, SARS-CoV-2)	EN 14476	30 sec.	60 sec.
BVDV virus	EN 14476	30 sec.	30 sec.
Rota virus	EN 14476	30 sec.	60 sec.
Noro virus	EN 14476	30 sec.	30 sec.
Hygienic hand disinfection	EN 1500	2 x 3 ml in 2 x 15 sec.	
Surgical hand disinfection	EN 12791	2 x 3 ml in 2 x 90 sec.	

Active ingredients: ethanol, propan-2-ol  
Packaging units: 500ml bottle, 1L bottle, 5L canister

Permission number for biocidal products in Poland: 7674/19

# Hand disinfection and hand care



### Logistic details

	500 ml	5L
Pcs / Box	24	1
No of boxes / pallet	40	96
Pcs / europallet	960	96

## Velodes® Gel

### Gel for hygienic and surgical hand disinfection

Velodes® Gel is a gel designed for hygienic and surgical hand disinfection. It has a bactericidal, yeastcidal, mycobacterial, and limited virucidal effect. It acts on enveloped viruses such as Vaccinia, BVDV, SARS-CoV-2, HIV, HBV, HCV, Rota and Noro viruses. Gentle for hands, it contains glycerine for skin care. Special formula prevents from hand stickiness. Dermatologically tested.

**Directions for use:** Hygienic hand disinfection (according to EN 1500): rub 2 x 3 ml of the product for 30 seconds (2 x 15 sec.). Surgical hand disinfection (according to EN 12791): rub 2 x 3 ml of the product for 3 minutes (2 x 90 sec.). Keep hands equally moist by adding more liquid.

	Gel		Active against B, F, V, Tbc, Rota, Noro		Hygienic hand disinfection in 30 sec.		Surgical hand disinfection in 3 min.
--	-----	--	---	--	---------------------------------------	--	--------------------------------------

Effectiveness	EN 14885	Conditions	
		clean	dirty
Bacteria (including MRSA)	EN 13727	30 sec.	30 sec.
Fungi (C. albicans)	EN 13624	15 sec.	15 sec.
Tuberculosis (M.terrae)	EN 14348	30 sec.	30 sec.
Enveloped viruses (Vaccinia , HIV, HBV, HCV, SARS-CoV-2)	EN 14476	30 sec.	60 sec.
BVDV virus	EN 14476	30 sec.	30 sec.
Rota virus	EN 14476	30 sec.	60 sec.
Noro virus	EN 14476	30 sec.	30 sec.
Hygienic hand disinfection	EN 1500	2 x 3 ml in 2 x 15 sec.	
Surgical hand disinfection	EN 14476	2 x 3 ml in 2 x 90 sec.	

Active ingredients: ethanol, propan-2-ol  
Packaging units: 500ml bottle, canister 5 L

Permission number for biocidal products in Poland: 7684/19

# Hand disinfection and hand care

## Logistic details

	500 ml	5 L
Pcs / Box	24	1
No of boxes / pallet	40	96
Pcs / europallet	960	96



pH 5,5 - 6,5

Contains: aqua, sodium laurethsulfate, cocamidopopyl betaine, glycerine, sodium chloride, lauramine oxide, panthenol, benzalkonium chloride, methylchloroisothiazolinone, methylisothiazolinone, citric acid

Packaging units: 500ml bottle, 5L canister

## Velodes® Scrub

### Antibacterial emulsion for hygienic hand washing

Velodes® Scrub is a delicate emulsion for washing hands and the whole body. It doesn't contain dyes or scented compositions. Perfect for doctors, dentists, nurses and people with sensitive skin, who wash their hands frequently. Neutral pH, the content of glycerine and coconut oil derivative allows everyday use without causing irritation. It contains panthenol, which gives anti-inflammatory effect and accelerates skin regeneration. It has properties unfavourable to bacterial growth. The emulsion can be used in the food industry, healthcare facilities, schools, kindergartens, and restaurants. Gives the skin sensation of softness. Dermatologically tested.

**Directions for use:** Apply a small amount of liquid on wet hands, rub to foam and rinse with water. Avoid contact with eyes



## Logistic details

	500 ml	5 L
Pcs / Box	24	1
No of boxes / pallet	40	96
Pcs / europallet	960	96



pH 5,5 - 6,5

Contains: aqua, glycerine, sodium chloride, sodium laurethsulfate, citric acid, cocamide dea, cocamidopropyl betaine, styrene/ acrylates copolymer, parfum, sodium lauryl sulfate, trideceth-7, methylchloroisothiazolinone, methylisothiazolinone, tetrasodium pyrophosphate, ascorbic acid

Packaging units: 500ml bottle, 5L canister

## Velodes® Soap

### Emulsion for hand and body washing

Velodes® Soap is an emulsion for washing hands before surgery and after hygienic hand disinfection. It can be also used for washing hands, hair and the whole body. The product does not contain soap, so it can be used by people with allergies and intolerance to soap-based products. Thanks to skin-nourishing ingredients, such as glycerine and coconut oil, it is particularly recommended for medical personnel susceptible to skin maceration related to frequent hand washing. The emulsion is also recommended for patients with bedsores. The product has a neutral pH. Dermatologically tested.

**Directions for use:** Hand washing: apply a small amount of the emulsion on wet hands and spread until it foams, then rinse with water. Washing in the shower: spread the emulsion on wet skin and hair. Wash your body, and rinse under running water. Washing in the bath: pour the emulsion into the bath filled with water. Wash your skin and hair, and then rinse under running water.



# Hand disinfection and hand care

## Logistic details

	1 L
Pcs / Box	6
No of boxes / pallet	80
Pcs / europallet	480



Contains: glycerine, sodium C14-16 olefin sulfonate, propylene glycol, cocamidopropylbetaine, lactic acid, methylchloroisothiazolinone, methylisothiazolinone  
 Packaging units: 1 L cartridge for dispenser

# Velodes® Soap Foam

## Foam for hand-washing

Velodes® Soap Foam is a gentle liquid applied in the form of foam for hygienic hand-washing. The product does not contain soap, so it can be used by people with allergies and intolerance to soap-based products. The liquid creates white, durable and thick foam which gently covers washed hands. Thanks to skin-nourishing and moisturizing ingredients, it is particularly recommended for medical personnel susceptible to skin maceration related to frequent hand washing. Dermatologically tested.

**Directions for use:** Hand washing: apply a small amount of the foam on wet hands, then rinse with water. Washing in the shower: spread the foam on wet skin and hair. Wash your body, and rinse under running water.

 Foam	 Hand washing
 Neutral scent	 White flowers scent



### Available models:



Elbow dispenser for closed system 1 L white  
 Product code: IM-DEZ-ML009  
 Dimensions: Width 133 mm, Height 270 mm, Depth 195 mm



Elbow dispenser for closed system 1 L white-violet.  
 Product code: IM-DEZ-ML009F  
 Dimensions: Width 133 mm, Height 270 mm, Depth 195 mm



Soap dispenser for closed system 1 L white.  
 Product code: IM-DEZ-ML007  
 Dimensions: Width 133 mm, Height 270 mm, Depth 115 mm



Soap dispenser for closed system 1 L white-violet.  
 Product code: IM-DEZ-ML007F  
 Dimensions: Width 133 mm, Height 270 mm, Depth 115 mm

# Wall-mounted - manual Dispensers

## Closed-system Dispenser

- Airtight closed-system with no air access
- Touch-free liquid application - prevents cross-infection
- Velodes 1 L disposable packaging in bag form
- Disposable pumps - no contamination risk
- Quick and hygienic cartridge refill
- Exceptionally easy to use and install

## ABS-plastic dispenser

- Dispenser in ABS plastic with high alcohol-resistance
- May be sterilised in an autoclave
- Has a movable disc allowing you to set the dose of 0.5 ml to 1.5 ml (every 0.5 ml)

### Available models:

- **ABS-plastic wall-mounted dispenser 500 ml.**  
Product code: IM-DEZ-ML011  
Dimensions: Width 85 mm, Height 280 mm, Depth 135 mm
- **ABS-plastic wall-mounted dispenser 1 L.**  
Product code: IM-DEZ-ML010  
Dimensions: Width 105 mm, Height 325 mm, Depth 160 mm
- **ABS-plastic wall-mounted dispenser 500 ml with stillage**  
Product code: IM-DEZ-ML012  
Dimensions: Width 158 mm, Height. 278 mm, Depth. 180 mm



## Dispenser in stainless steel

The stainless-steel elbow dispenser is a combination of minimalism, simplicity, and aesthetics, with the robustness of the materials used.

### Available models:

- **Stainless-steel wall-mounted dispenser 500 ml.**  
Product code: IM-DEZ-ML020D  
Dimensions: Width 77 mm, Height 198 mm, Depth 55 mm
- **Stainless-steel wall-mounted dispenser 1 L.**  
Product code: IM-DEZ-ML021D  
Dimensions: Width 95 mm, Height 225 mm, Depth 65 mm



## Dispenser stands PROCLEAN ECO

Stands with H and O bases allowing PROCLEAN 05 ECO and 1 ECO stainless steel dispensers to be mounted anywhere.

### Available models:

- **Stand with H base**  
Product code: IM-DEZML028D  
Dimensions: width 360 mm, height 1115,5 mm, depth 310 mm
- **Stand with O base**  
Product code: IM-DEZ-ML031D  
Dimensions: width 1115,5 mm, diameter 350 mm



Stands fitted for dispensers in stainless steel

## Automatic, stand-alone and wall-mounted dispensers

- Touch-free operation- the device automatically activates if a hand is brought closer to the dispenser
- Infrared motion detector - the special motion detector is also suitable for use in poor lighting conditions
- Automatic liquid dispensing - the device measures the amount of disinfectant to be dispensed
- Robust metal design - disinfectant is in impact-resistant metal
- Mains and battery power supply (2in1) - the dispenser can be connected to a 230V socket or powered by 8 x AA batteries
- Liquid container - the dispenser design is compatible with the replaceable cartridges of Velodes liquids

### Available models:

- Product code: IM-DEZ-ML024D  
Dimensions: width 280 mm, height 1301 mm, depth 240 mm
- Product code: IM-DEZ-ML033D  
Dimensions: width 155 mm, height 323 mm, depth 85 mm



IM-DEZ-ML024D

IM-DEZ-ML033D

# Hand disinfection and hand care



### Logistic details

	100 ml	500 ml
Pcs / Box	50	30
No of boxes / pallet	72	40
Pcs / europallet	3600	1200

Ingredients: aqua, cetostearylalcohol, glycerine, glyceryl stearate, isopropyl palmitate, cetareth-20, cyclopentasiloxane, paraffinumliquidum, phenoxyethanol (and) ethylhexylglycerine, DMDM hydantoin, carbomer, parfum, citric acid, BHT, oleaeuropaea (olive) oil, tocopherylacetate, triethanolamine, linalool, citronellol, alpha-isomethylionone, geraniol, limonene, hexyl cinnamal

Packaging units: 100 ml tube, 500 ml bottle

## Velodes® Cream

### Protective hand cream

Velodes® Cream is a gentle cream for hands and body skin care. Protects, moisturizes and regenerates. It nourishes dry, irritated and rough skin. It is recommended for everyday use after frequent hand washing and hygienic disinfection. It is particularly recommended for medical personnel and people exposed to drying and maceration of the skin, as a result of frequent washing and wearing gloves. Provides silky hands softness. Quick absorption brings an immediate feeling of comfort. It does not contain parabens. Dermatologically tested.

**Directions for use:** Apply a small amount of cream on clean, dry hands. Gently rub into skin of your hands and nails. Repeat the process after each hand washing.

Cream	Moisturizes, protects and cares for hands
Contains olive oil	Contains vitamin E and glycerine



### Logistic details

	100 ml	500 ml
Pcs / Box	50	30
No of boxes / pallet	72	40
Pcs / europallet	3600	1200

Ingredients: aqua, cetostearylalcohol, glycerine, glyceryl stearate SE, isopropyl palmitate, cetareth-20, phenoxyethanol, aloe barbadensisleaf extract, cyclopentasiloxane, paraffinumliquidum, D-panthenol, carbomer, parfum, BHT, tocopherylacetate, triethanolamine, ethylhexylglycerine, benzyl salicylate, limonene, butylphenyl methylpropional

Packaging units: 100 ml tube, 500 ml bottle

## Velodes® Cream Aloe Vera

### Protective hand cream


Velodes® Cream Aloe Vera moisturizes, nourishes the skin of hands and protects it from the adverse influence of external factors. Rich formula, containing glycerine, aloe extract and vitamin E, makes hand skin soft and silky. It is recommended for everyday use after frequent hand washing and after hygienic disinfection. It is particularly recommended for medical personnel and people exposed to drying and maceration of the skin, as a result of frequent washing and wearing gloves. Dermatologically tested.

**Directions for use:** Apply a small amount of cream on clean, dry hands. Gently rub into skin of your hands and nails. Repeat the process after each hand washing.

Cream	Moisturizes, protects and cares for hands
Contains aloe extract	Contains vitamin E and glycerine

# Hand disinfection technique

Efficient protection against bacteria, mycobacteria, yeasts and viruses: HIV, HBV, HCV, Ebola, influenza. Extended activity time: up to 3 hours.




Take an appropriate amount of the antiseptic product Velodes® Soft.

Rub one palm against the other palm. Each stage involves five movements in both directions.

Rub the dorsal part of one hand with the palm of the other hand and change hands.

Rub palm to palm with interlaced spread fingers.



Rub the backs of closed fingers of one hand against the bent fingers of the other hand and change hands.

Take your right thumb with your left hand and perform rotational movements around the thumb and change hands.

Rub the fingertips of the right hand on the left palm with circular movements and change hands.

After drying hands are disinfected.



# MEDISEPT

back



## MEDISEPT

We are a leading European manufacturer of disinfectants, and professional cleaning agents. Our products are used in and by hospitals, private clinics, dental and medical surgeries, public institutions, and are available in medical shops and via network of distributors.

Medisept has implemented a Quality Management System for medical manufacture and has been certified in accordance with EN ISO 13485 in the scope of design, production and distribution of invasive and non-invasive medical products.

In 2014, we launched the most modern production line in the industry and an innovative research centre in factory.

We effectively protect people by providing modern and comprehensive solutions in the field of hygiene and disinfection – this is our mission, that we are implementing in our everyday work. Our values help us to achieve this goal: commitment, cooperation and knowledge sharing, responsibility, professionalism, consistency in action, and respect.



Manufacturer:  
MEDISEPT Sp. z o.o.,  
Konopnica 159 c  
21-030 Motycz, Poland

tel.: +48 81 535 22 66  
fax: +48 81 535 22 37  
export@medisept.pl

rev. 2021.07

[www.medisept.pl](http://www.medisept.pl)