



MEDISEPT

Perfect disinfection

Perfect disinfection

Catalogue 2023



MEDISEPT - European company for You!

About as

MEDISEPT is a European producer of comprehensive solutions in the field of hygiene and disinfection in the medical standard. We have been operating on the market since 1994, providing solutions for the professional sector, including health care, industry and cleaning companies.

Why is it worth to choose MEDISEPT products?

All our products are a response to the requirements of the medical market and meet the highest world quality standards. Biocides offered by our company are manufactured to the same quality standards as medical disinfectant products. MEDISEPT products guarantee high efficiency and an unique attractive price alternative while maintaining incomparable quality.

The offer of MEDISEPT

We offer the highest quality disinfectants which for years have been able to effectively fight pathogenic microorganisms. They can also help to maintain personal safety both in the hospital environment and in public institutions. All MEDISEPT products are a complementary offer developed in accordance with the system of norms and certificates regulating the standards on the professional market.

Standards and certification

Quality Management System for medical devices

All MEDISEPT disinfection products and cleaning agents are manufactured in accordance with the ISO 13485¹ quality assurance system for medical devices certified by TUV Nord². Each product is tested according to the latest harmonized standards confirming the microbiological effectiveness against the declared microorganisms³.



1 EN ISO 13485 is the international standard that defines the requirements for quality management systems, the fulfillment of which is a decisive factor that builds the trust of customers in the medical sector. Thanks to the constant monitoring of all key processes in the company, we guarantee the highest standards of production and service to our customers not only in the medical industry, but also in other market segments.

2 Meeting legal requirements means that a medical device is designed, manufactured and monitored at every stage of the product life cycle. The release of each lot for sale is preceded by the control of quality parameters, and its fulfillment is confirmed by a lot analysis certificate. Placing on the market is carried out on the basis of suitable marketing authorizations and registration in appropriate offices of a given country.

3 The effectiveness of disinfectants is confirmed by a set of tests defined in the EN 14885 standard, carried out only in accredited research laboratories.

Discover MEDISEPT disinfection products

We present a line of preparations for disinfecting surfaces, tools and hands, intended for both medical and non-medical use. Disinfectants remove harmful microorganisms, bacteria, viruses or tuberculosis mycobacteria. Preparations not only eliminate the risk, but also help to protect the surface and tools. Hand disinfectants, in addition to providing protection against microorganisms, affect their regeneration processes. We offer ready-to-use liquids as well as powders or concentrates for disinfection. MEDISEPT's offer includes both products for everyday cleaning of surfaces and tools, as well as alcohol-free sprays and wipes that allow the neutralization of microorganisms in any place.



How to read labels?

All MEDISEPT product's labels contain four sections of pictograms always arranged in the same order. Concentrates have the additional field with dosing instructions. Each section is designed to visually convey the most important information about the product.



1. Section: product form

Pictograms informing about the form of the product. It can be spray, foam, powder, liquid, emulsion, concentrate or wipes.

2. Section: distinguishing feature

Pictograms indicating the use of the product, disinfection, quick disinfection, washing or care.

3. Section: destiny

Pictograms indicating the type of surfaces the product is intended for.



Neutral
scent

Dosage:

0,5% → 5 ml in 995 ml of water 15 min.
1% → 10 ml in 990 ml of water 15 min.

4. Section: fragrance

Pictograms and description on the background of the color corresponding to the smell of the product.

5. Section: dosing

This section provides dosing guidance for each concentration and time.

The main features of the packaging

Integrated dispenser
for the easy dosing of
a concentrate

High-quality white
plastic

The background
colour of the label
corresponding to the
application

Product name

Brief description of
the product

Sections for product-
form Pictograms,
target areas and
additional
information

Dosage instructions

Product-category
code



information about
efficacy

Information about
dermatological tests

Marketing
authorization number

Tolerance to sensitive
materials

QR code which
redirects to the product
website

Marking medical device

Colour-coding of products

Target areas:

SURFACES

INSTRUMENTS

HANDS

Intended use category:

DISINFECTION

CARE / WASHING

10 benefits, advantages and profits of MEDISEPT disinfection products

Marketing authorization issued by the relevant institutions



1

Full microbiological effectiveness (viruses, bacteria, fungi, mycobacteria) performed in accredited laboratories



2

Dermatologically ingredients safe for skin in long term contact and pleasant fragrance



3

Experience in creating safe formulations



4

Production quality control according to ISO standard



5

Experience in developing procedures, hygiene plans for all industries



6

Material tolerance confirmed by utility tests



7

Full range of disinfectant products: for hands, surfaces and instruments in various forms of use



8

Production of biocidal products and cleaning agents in medical standard



9

Many years of experience in working with hospitals and isolation wards - Knowledge of medical procedures



10





Additional information



Broad microbicidal spectrum



sporocidal

NEW

Logistic details

Cap.	Product code	Scent	Pcs / Box	No of boxes / pallet	Pcs / europallet
1L	SSE-ALKGOT-ML779	-	12	8	384
5L	SSE-ALKGOT-ML526	-	1	30	120

Velox® Oxy ETA

Spray for disinfection and surface cleaning

Product features:

- ▶ Broad sporocidal activity in 15 min. (*C. difficile*, *B. subtilis*)
- ▶ Full fungicidal activity in 30 sec. (*A. brasiliensis*, *C. albicans*)
- ▶ Active against enveloped viruses (Vaccina, HIV, HBV, HCV, SARS-CoV-2) in 30 sec and non-enveloped viruses (Noro) in 15 sec.

Directions for use:

Spray the surface from the distance of about 20-30 cm, ensuring it is fully covered. Remove the excess of the preparation using e.g. a disposable towel. Surface does not require rinsing with water. Before disinfecting particularly sensitive surfaces (e.g. natural leather), it is recommended to perform the test on a non-visible part of the surface. For professional use only. Do not mix with other cleaning products or disinfectants. Read the information leaflet before use.

Application:

Velox® Oxy ETA is recommended for the disinfection of small surfaces of medical equipment, including treatment beds and chairs, patient cabinets, medical and operating equipment, incubators, rehabilitation equipment, as well as dental tips, units, treatment lamps, worktops, prosthetic trays, prostheses and impressions. Velox® Oxy ETA can be used on alcohol resistant materials.



Spray



Active against Polio, Adeno, Noro, Spores



Medical and non-medical surfaces



Hairdressing and beauty tools

Active substances in 100g of product:

72,8 g Ethanol, 7,2 g, Propan-2-ol, 3 g hydrogen peroxide

Effectiveness

Bacteria

Fungi (*C. albicans*)

Fungi (*A. brasiliensis*)

Mycobacterium tuberculosis (*M.terrae*, *M.avium*)

Enveloped viruses

Adeno Virus

Noro Virus

Polio Virus

Spores (*B. subtilis*)

Spores (*Cl. difficile* R027)

EN 14885

Conditions

clean

dirty

EN 13727

EN 16615

EN 13697

EN 13624

EN 16615

EN 13697

EN 13624

EN 13697

EN 16615

EN 14348

EN 14476

EN 14476

EN 14476

EN 14476

EN 13704

EN 17126

60 sec.

15 sec.

30 sec.

15 sec.

30 sec.

60 sec.

30 sec.

60 sec.

30 sec.

3 min

30 sec.

30 sec.

15 sec.

15 sec.

10 min

15 min



Recommended products

- ▶ MEDISEPT Dry Wipes → p.15
- ▶ MEDISEPT Dry Wipes Box → p.15



Matching products

- ▶ IM-AKC-0033 → p.34

Packaging units: 1L spray bottle, 5L canister

CE 2274 class II a



Available scents:



Neutral
scent



Grapefruit
scent

Logistic details

Cap.	Product code	Scent	Pcs / Box	No of boxes / pallet	Pcs / europallet
1L	SSE-43-ALKGOT-ML082	neutral	12	8	384
	SSE-43-ALKGOT-ML082N	neutral			
	SSE-43-ALKGOT-ML924	grapefruit			
	SSE-43-ALKGOT-ML924N	grapefruit			
5L	SSE-43-ALKGOT-ML099	neutral	1	30	120
	SSE-43-ALKGOT-ML931	grapefruit			

N - bottle with sprayer
without N - bottle without sprayer

Velox® Top AF

Spray for disinfection and surface cleaning

Product features:

- ▶ Active in 30 sec. against Bacteria, Fungi, Tuberculosis, enveloped viruses, Adeno, Rota
- ▶ Recommended for small surfaces of medical equipment, chairs, beds, or other difficult to access surfaces, as well as food related objects. Also recommended for disinfection of dental tips before sterilization and of prosthesis and silicone impressions
- ▶ Intended for all non-medical surfaces e.g. gastronomy, hairdressing salons and all surfaces having contact with food
- ▶ Free from aldehyde and phenol, does not discolour disinfect
- ▶ High material tolerance to ABS and upholstery materials - Famed opinion
- ▶ Available in two scent versions
- ▶ Possible registrations: medical device and biocide

Directions for use:

Spray the surface from the distance of 30 cm covering the whole disinfected area. Spread the liquid using disposable towels and allow to dry. The surface does not need rinsing with water, dries out completely and does not leave smudges. Do not use on alcohol sensitive surfaces (i.e. acrylic glass). For professional use only.

Application:

Velox® Top AF is an alcohol-based liquid for cleaning and rapid disinfection of small surfaces of treatment chairs, operating tables, medical equipment and dental handpieces before sterilization process. Also intended for use in gastronomy, nursing homes, sanatoriums, fitness and SPA salons, salons hairdressers, beauty salons, massage rooms, kindergartens, schools and for disinfecting surfaces having contact with food. It is ideal for cleaning and disinfection of small and hard-to-reach surfaces such as countertops, tables, handles, handles, and for hairdressing and cosmetic equipment and accessories.



Spray



Active against
Tbc, Rota,
Adeno



Active in
just 30 sec.



Medical and
non-medical
surfaces

Active substances in 100g of product:

40 g propan-2-ol, 20 g ethanol, 0.39 g amines, n-C10-16-alkyltrimethyl-di-products reaction with chloroacetic acid, 0.25 g didecyldimethylammonium chloride

Effectiveness

	EN 14885	Conditions	
		clean	dirty
Bacteria	EN 1276	30 sec.	60 sec.
	EN 13727	30 sec.	60 sec.
	EN 13697	30 sec.	60 sec.
Fungi (C. albicans)	EN 1650	30 sec.	60 sec.
	EN 13624	30 sec.	60 sec.
Mycobacterium tuberculosis (M. terrae)	EN 14348	30 sec.	60 sec.
Enveloped viruses	DVV/RKI	30 sec.	30 sec.
Rota Virus	EN 14476	30 sec.	
Adeno Virus	EN 14476	60 sec.	
BVDV Virus	DVV/RKI	30 sec.	



Neutral



Grapefruit



Recommended products

- ▶ MEDISEPT Dry Wipes → p.15
- ▶ MEDISEPT Dry Wipes Box → p.15



Matching products

- ▶ IM-AKC-0033 → p.34

Packaging units: 1L spray bottle, 5L canister

CE 2274 class IIa



Available scents:



Neutral
scent



Tea Tonic
scent

Neutral scent



Flächendesinfektion ohne Mechanik / Surface disinfection without mechanical action

Wischen / Mechanical action		Organische Belastung / Interfering substance		Einwirkzeit in min / Contact time in min			
mit / with	ohne / without	gering / clean conditions	hoch / dirty conditions	5	15	30	60 240
	x		x	konz.			

Logistic details

Cap.	Product code	Scent	Pcs / Box	No of boxes / pallet	Pcs europallet
1L	SSE-ALKGOT-ML680	neutral	12	8	384
	SSE-ALKGOT-ML680N	neutral			
	SSE-ALKGOT-ML703	tea tonic			
	SSE-ALKGOT-ML703N	tea tonic			
5L	SSE-ALKGOT-ML697	neutral	1	30	120
	SSE-ALKGOT-ML710	tea tonic			

N - bottle with sprayer
without N - bottle without sprayer

Velox Spray

Spray for disinfection and surface cleaning

Product features:

- For disinfection of small surfaces of medical equipment: beds and treatment chairs, medical and surgical equipment, rehabilitation equipment
- Designed for all non-medical surfaces e.g. in gastronomy, hairdressing and for surfaces having contact with food
- Active already in 15 sec. against Fungi and in 30 sec. against Bacteria, Tuberculosis and viruses (including Rota, Noro)
- Does not leave streaks nor stains
- Available in two scent versions
- Possible registrations: medical device and biocide

Directions for use:

Spray the surface from the distance of about 30 cm, ensuring it is fully covered. Leave it for 30 seconds. In case of disinfecting surfaces having contact with food, rinse them with drinking water. Other surfaces do not require rinsing with water. Do not use on alcohol sensitive surfaces (eg. acrylic glass).

Application:

Alcoholic liquid for rapid disinfection of small areas in the medical area, nursing homes, sanatoriums, fitness and SPA centres, hairdressers, massage and beauty salons, kindergartens and schools. The product is suitable for surfaces having and not-having contact with food. It can be used for disinfecting surfaces, such as table tops, tables, handles, for hairdressing and cosmetics equipment and accessories. The product is recommended also for disinfection of non-invasive medical devices such as treatment chairs and medical equipment. Velox® Spray is also recommended for the disinfection of dental handpieces (before sterilisation), prosthetic trays and silicone dental impressions.



Spray



Surfaces
disinfection



Hairdressing
and cosmetic
tools



Non-medical
surfaces

Active substances in 100g of product:

63,7 g ethanol, 6,3 g propan-2-ol

Effectiveness

EN 14885

Conditions

clean

dirty

Bacteria

EN 13727

30 sec.

30 sec.

EN 16615

30 sec.

30 sec.

EN 13697

30 sec.

60 sec.

Fungi (C. albicans)

EN 13624

15 sec.

15 sec.

EN 16615

30 sec.

30 sec.

EN 13697

30 sec.

Mycobacterium tuberculosis (M. terrae)

EN 14348

30 sec.

30 sec.

Enveloped viruses

EN 14476

30 sec.

60 sec.

BVDV Virus

EN 14476

30 sec.

30 sec.

Rota Virus

EN 14476

30 sec.

60 sec.

Noro Virus

EN 14476

30 sec.

30 sec.

Bactericidal activity and yeasticidal activity

VAH

5 min



Neutral



Tea tonic



Recommended products

- › MEDISEPT Dry Wipes → p.15
- › MEDISEPT Dry Wipes Box → p.15



Matching products

- › IM-AKC-0033 → p.34

Packaging units: 1L spray bottle, 5L canister

CE 2274 class IIa

Check
availability

Logistic details

Cap.	Product code	Scent	Pcs / Box	No of boxes / pallet	Pcs / europallet
Box	IM-DEZ-ML979	-	6 or 12	10 or 19	912 or 960
Refill	IM-DEZ-ML962	-	12	16 or 24	1344 or 2304

Velox® Wipes

Alcohol-based wipes for cleaning and disinfection of small surfaces and medical equipment

Product features:

- ▶ Active in 30 secs against Bacteria, Fungi, Tuberculosis, enveloped viruses, Adeno and Rota viruses
- ▶ Recommended for units, medical equipment, rehabilitation equipment, dental handpieces, treatment lamps, patient beds
- ▶ Suitable for medical equipmentsurfaces having contact with food
- ▶ High material tolerance to ABS and upholstery materials - Famed opinion
- ▶ Non-woven wipes 13x20 cm, 23 g/m²
- ▶ Possible registrations: medical device and biocide

Directions for use:

Place the refill in the box. Take a tissue out of the box and wipe a surface to be cleaned and disinfected. Close the lid. After opening the package, use the tissues within 21 days. Do not use if a tissue is dry. Do not use on areas sensitive to alcohol (e.g. acrylic glass). For professional use only.

Application:

Velox® Wipes are dedicated for rapid disinfection of small, difficult to access surfaces of medical equipment, treatment chairs, beds, dental ends before the sterilisation process.



Wipes

Active against
Adeno, RotaMedical and
non-medical
surfacesActive
in 30 sec.

Active substances in 100g of soaking solution

40 g propan-2-ol, 20 g ethanol, 0.39 g amine, n-C10-16 alkyltrimethylenedio di-reaction products with chloroacetic acid, 0.25 g didecyl-dimethylammonium chloride

Effectiveness

	EN 14885	Conditions	
		clean	dirty
Bacteria	EN 1276	30 sec.	60 sec.
	EN 13697	30 sec.	60 sec.
	EN 13727	30 sec.	60 sec.
	EN 16615	30 sec.	30 sec.
Fungi (C. albicans)	EN 1650	30 sec.	60 sec.
	EN 13624	30 sec.	60 sec.
	EN 16615	30 sec.	30 sec.
Mycobacterium tuberculosis (M. terrae)	EN 14348	30 sec.	60 sec.
Enveloped viruses	DVV/RKI	30 sec.	30 sec.
Rota Virus	EN 14476	30 sec.	
Adeno Virus	EN 14476	60 sec.	
BVDV Virus	DVV/RKI	30 sec.	



Packaging units: 100 pieces refill., 100 pieces box

Size: 13 x 20 cm

Grammage: 23 g/m²

CE 2274 class IIa

Available scents:



Tea Tonic
scent



Logistic details

Cap.	Product code	Scent	Pcs / Box	No of boxes / pallet	Pcs / europallet
Box	IM-DEZ-ML588	Tea tonic	6 or 12	10 or 19	912 or 960
Refill	IM-DEZ-ML595		12	16	1344
Flow pack	IM-DEZ-ML410		15	9	945

Velox® Duo Wipes

Alcoholic wipes for cleaning and rapid disinfection of small surfaces and medical equipment

Product features:

- ▶ Active in 30 secs against Bacteria, Fungi, Tuberculosis, enveloped viruses, Adeno and Rota viruses
- ▶ Recommended for medical equipment, rehabilitation equipment, dental handpieces, patient beds, cabinets, tabletops
- ▶ For disinfection of medical devices surfaces having contact with food
- ▶ Non-woven wipes 19x15 cm, 50 g/m² - flow pack
- ▶ Non-woven wipes 13x20 cm, 23 g/m² - box / refill
- ▶ Possible registrations: medical device and biocide
- ▶ Positive opinion of Famed

Directions for use:

Take a wipe out of the package and wipe the surface. Wait 30 seconds. In dirty conditions operating time against enveloped viruses and Rota virus is extended up to 60 seconds. In case disinfecting surfaces having contact with food, rinse them with potable water. Other surfaces do not require rinsing with water. Close the package tightly. After opening the package, use the wipes within 21 days. Do not use if the wipe is dry. Do not use on surfaces sensitive to alcohol (e.g. acrylic glass). Product intended for use by the general user.

Application:

Wipes for cleaning and rapid disinfection of small surfaces in medical field, nursing homes, sanatoriums, fitness and SPA centres, hairdressing, massage and beauty salons, kindergartens and school. The product is suitable for surfaces having contact with food. Recommended for small surfaces of medical equipment: operating beds and chairs, beds, medical and surgical equipment and disinfection of dental tips before sterilisation process, prosthetic spoons and silicon moulds.



Wipes



Active in
30 sec.



Medical and
non-medical
surfaces



Surfaces
having contact
with food

Active substances in 100g of soaking
solution

63,7 g ethanol, 6,3 g propan-2-ol

Effectiveness	EN 14885	Conditions	
		clean	dirty
Bacteria	EN 13727	30 sec.	30 sec.
	EN 16615	30 sec.	30 sec.
	EN 13697	30 sec.	60 sec.
Fungi (C. albicans)	EN 13624	15 sec.	15 sec.
	EN 16615	30 sec.	30 sec.
	EN 13697	30 sec.	
Mycobacterium tuberculosis (M. terrae)	EN 14348	30 sec.	30 sec.
Enveloped viruses	EN 14476	30 sec.	60 sec.
BVDV Virus	EN 14476	30 sec.	30 sec.
Rota Virus	EN 14476	30 sec.	60 sec.
Noro Virus	EN 14476	30 sec.	30 sec.



Packaging units: 50 pcs, flow pack with adhesive label, Size: 19 x 15 cm, Grammage: 50 g/m²
Packaging units: Box 100 pcs, refill 100 pcs, Size: 13 x 20 cm, Grammage: 23 g/m²

CE 2274 class IIa



Available scents:



Cashmere
scent

Logistic details

Cap.	Product code	Scent	Pcs / Box	No of boxes / pallet	Pcs europallet
Flow pack	IM-DEZ-ML976	Cashmere	15	9	945

MEDISEPT Wipes for hand and surface disinfection

High-quality wipes for disinfecting hands and all small, alcohol- resistant surfaces

Product features:

- ▶ Tested in accordance with new standards for disinfectants
- ▶ Hygienic disinfection in 30 sec. (soaking liquid tested in acc. with EN 1500)
- ▶ High quality, non-woven fabric: 50 g/m²
- ▶ The capacity of one wipe is 0.5 m with a double rubbing*
- ▶ Non-woven: 19 x 15 cm
- ▶ Possible registrations: medical device and biocide

Directions for use:

Take a wipe out of the package and wipe hands or surface. Close the package tightly. After opening, use the wipes within 21 days. Do not use if the wipe is dry.

Application:

MEDISEPT Wipes for hand and surface disinfection are alcoholic wipes designed for hygienic hands disinfection and any small alcohol-resistant surfaces. They are recommended for surfaces such as tabletops, knobs, handles, seating surfaces, medical equipment, and accessories made of plastic, stainless steel, porcelain and ABS plastic.



Wipes



Works in 30 sec.
on Tbc. and Viruses



Disinfect
and protects
hands



Disinfect
surfaces

Active substances in 100g of soaking solution

60 g propan-2-ol, 0,5 g glukonian chlorheksydyny

Effectiveness

	EN 14885	Conditions	
		clean	dirty
Bacteria	EN 13727	30 sec.	60 sec.
	EN 16615		30 sec.
Fungi (C. albicans)	EN 1275	30 sec.	
	EN 16615		30 sec.
	EN 13624	30 sec.	
Mycobacterium tuberculosis (M. terrae)	EN 14348	30 sec.	
Enveloped viruses	EN 14476	30 sec.	
Rota Virus	EN 14476	45 sec.	
BVDV Virus	EN 14476	30 sec.	
Hygienic hands disinfection*	EN 1500	2 x 1.5 ml in 2 x 15 seconds	

*Test performed on saturated formulation, 1 tissue contains approx. 3.56 ml of disinfecting liquid



Packaging units: 50 pcs. flow pack with self-adhesive closing

Size: 19 x 15 cm

Grammage: 50 g/m²

CE 2274 class IIa

Disinfection of surfaces without chlorine and alcohol



Non-alcoholic preparations, based on active substances such as QAV, amines have a microbicidal effect against enveloped viruses, which include SARS-CoV-2 virus.

Benefits of non-alcoholic products



Efficiency and economy

Specially designed for large area disinfection - 5L of concentrate is sufficient for up to 20.000 m²

Safety

They do not cause irritating fumes - i. e. in the case of alcohol and chlorine - therefore can be used over large areas in the presence of people. They act on the full spectrum: viruses, bacteria, fungi, mycobacteria

Do not damage the surface

Suitable for disinfecting alcohol-sensitive surfaces such as wood, acrylic, leather, plastics. They are approved for material tolerances for medical devices



Methods of application of non-alcoholic disinfection



Spraying of the surface with the sprayer



Floor cleaning with a mop



Wiping the surface with a cloth



Logistic details

Cap.	Product code	Scent	Pcs / Box	No of boxes / pallet	Pcs / europallet
1L	SSE-GOT-ML405	-	12	8	384
5L	SSE-GOT-ML412	-	1	30	120

Velox® Foam Extra

Foam for cleaning and quick disinfection of sensitive surfaces of non-invasive medical devices

Product features:

- ▶ Alcohol-free foam
- ▶ Active in 1 min. against Bacteria, Fungi, Viruses and TBC in 5 min
- ▶ Ready to use, for surfaces sensitive to alcohol
- ▶ Recommended for disinfection of various surfaces of medical equipment made of glass, porcelain, metal, rubber, plastic and acrylic glass. Also for treatment chairs, rehabilitation equipment and incubators
- ▶ High material tolerance for ABS and upholstery materials - Famed opinion
- ▶ Pleasant scent
- ▶ Possible registrations: medical device and biocide

Directions for use:

The product has a double, foam and spray head. Application in foam is recommended. Spray the surface from a distance of 20-30 cm, making sure it is fully covered. Remove the excess foam with a disposable paper towel. The surface does not need to be rinsed with water. Before disinfecting particularly sensitive surfaces (e.g. natural leather), it is recommended to perform the test on a non-visible part of the surface. For professional use only.

Application:

Foam for disinfecting a wide range of surfaces in medical areas, nursing homes, fitness and spa salons, beauty and massage parlors. It can be used in intensive care unit and operating theater, for disinfection of medical equipment, rehabilitation equipment, treatment chairs, incubators and treatment lamps. Does not contain aldehydes and phosphates. Recommended for alcohol-sensitive surfaces such as acrylic glass. It has pleasant scent. Can be used for surfaces in contact with food. The product has a positive opinion of medical equipment manufacturer Famed.



Foam



Active against Tbc



Medical and non-medical surfaces



Alcohol-free product

Active substances in 100g of product:

0,15 g N-(3-aminopropyl)-N-dodecylpropano-1,3-diamina, 0,14 g poli(oksy-1,2-etanodilo), alfa.-[2-(didecylmetyloamino)etylo]-.omega.-hidroksy-.propanian(sól)

Effectiveness

	EN 14885	Conditions	
		clean	dirty
Bacteria	EN 13727	60 sec.	60 sec.
	EN 16615	60 sec.	3 min
	EN 13697	30 sec.	60 sec.
Fungi (C. albicans)	EN 13624	60 sec.	60 sec.
	EN 16615	60 sec.	60 sec.
	EN 13697	30 sec.	
Mycobacterium tuberculosis (M.terrae)	EN 14348	5 min	10 min
Mycobacterium tuberculosis (M. avium, M. terrae)	EN 14348	15 min	15 min
Enveloped viruses	EN 14476	60 sec.	60 sec.
BVDV Virus	EN 14476	60 sec.	60 sec.



Recommended products

- ▶ MEDISEPT Dry Wipes → p.15
- ▶ MEDISEPT Dry Wipes Box → p.15



Matching products

- ▶ IM-AKC-0033 → p.34

Packaging units: 1L spray bottle, 5L canister

CE 2274 class IIa



Logistic details

Cap.	Product code	Scent	Pcs / Box	No of boxes / pallet	Pcs / europallet
Box	IM-DEZ-ML993	-	6 or 12	10 or 19	912 or 960
Refill	IM-DEZ-ML986	-	12	16 or 24	1344 or 2304

Velox® Wipes NA

Alcohol-free wipes for cleaning and rapid disinfection of small surfaces and medical equipment

Product features:

- ▶ Active in 1 min against Bacteria, Fungi, enveloped viruses
- ▶ Active in 5 min against Tuberculosis
- ▶ Recommended for ultrasound tubes and heads
- ▶ Ready to use, for surfaces sensitive to alcohol
- ▶ Possible registrations: medical device and biocide
- ▶ Pleasant scent
- ▶ Destined for disinfection of various medical equipment surfaces made of glass, porcelain, steel, rubber, plastic and acrylic glass, operating chairs, rehabilitation equipment, incubators
- ▶ Non-woven wipes 13 x 20 cm, 23 g/m²
- ▶ Positive opinion of famed

Directions for use:

Place the refill in the box. Take a tissue out of the box and wipe a surface to be cleaned and disinfected. Close the lid. After opening the package, use the tissues within 30 days. Do not use if a tissue is dry. The surface does not need to be rinsed with water. Before applying on particularly sensitive surfaces, test the wipes on a non-visible area. For professional use only.

Application:

Velox® Wipes NA are recommended for disinfecting surfaces of noninvasive medical equipment made of plastic, acrylic glass, stainless steel, metal, aluminium, rubber and porcelain. Suitable for disinfection of medical appliances, rehabilitation equipment, treatment chairs, incubators, the heads of ultrasound machines and lamps in intensive care units and operating theatres. Free from alcohol, aldehydes and phosphates, the wipes do not discolour the surfaces.



Wipes



Active in 1 min.
against B, F, V



Medical and
non-medical
surfaces



Alcohol-free
product

Active substances in 100g of soaking solution

0.15 g N- (3-aminopropyl)-N-dodecylpropane-1,3-diamine, 0.14 g Poli (oxy-1,2-ethanedyl), alpha - [2-(didecylmethylamino) ethyl] - omega. hydroxy, propanoate (salt)

Effectiveness

	EN 14885	Conditions	
		clean	dirty
Bacteria	EN 13727	60 sec.	60 sec.
	EN 16615	60 sec.	3 min
	EN 13697	30 sec.	60 sec.
Fungi (C. albicans)	EN 13624	60 sec.	60 sec.
	EN 16615	60 sec.	60 sec.
	EN 13697	30 sec.	
Mycobacterium tuberculosis (M. terrae)	EN 14348	5 min	10 min
Mycobacterium tuberculosis (M. avium, M. terrae)	EN 14348	15 min	15 min
Enveloped viruses	EN 14476	60 sec.	60 sec.
BVDV Virus	EN 14476	60 sec.	60 sec.



Packaging units: 100 pieces refill., 100 pieces box

Size: 13 x 20 cm

Grammage: 23 g/m²

CE 2274 class IIa



Logistic details

Cap.	Product code	Scent	Pcs / Box	No of boxes / pallet	Pcs / europallet
Bucket with wipes	IM-DEZ-ML035	-	1	24	24
Refill	IM-DEZ-ML059	-	24	4	672

MEDISEPT Dry Wipes / MEDISEPT Dry Wipes Box

Dry wipes for soaking with disinfecting / cleaning agent

Product features:

- ▶ The porous structure of the wipes with a grammage of 45 g/m² perfectly absorbs the agent and guarantees its ideal distribution on the cleaned surface
- ▶ The thick unwoven cloth does not scratch the surface
- ▶ Non-dusting and lint-free white wipes
- ▶ Very enduring and resistant to detergents and solvents
- ▶ Recommended soaking: 1 litre per bucket allows for effective 2 disinfection of 2 m surface with one wipe
- ▶ The wipe size is 16 x 30 cm.
- ▶ The optimal number of wipes (100 pcs) allows their effective use for up to 30 days (solution-storage time)
- ▶ Box with refill

Directions for use:

Prepare the soaking liquid according to the recommendations on the product's label. To soak the wipes in the liquid, follow the instructions as shown in the pictograms below. Wipe the surface with a damp wipe and leave it to dry for the time recommended by the soaking-liquid's manufacturer. Clean the bucket and lid before re-use with the soaking solution. Using the liquid-soaked wipes does not require protective gloves

1

Open the cover of the bucket. Unwrap the roll from the foil, throw away the foil.

2

Remove the first cloth from the center of the roll. Put it through the hole in the cover.

3

Fill the bucket with the selected solution of the product for cleaning / Disinfection with the desired concentration or product ready for use. Amount of soaking liquid: 1 liter.

4

Close the bucket and the opening. Leave for 30 minutes so that all the cloths have been fully soaked.

5

Fill in the sticker on the bucket. If you are using a refill, put a sticker on the bucket.

Application:

MEDISEPT Dry Wipes are designed for soaking with a ready-to-use product or utility solution of concentrate with water. Formulated for cleaning and disinfecting various surfaces of non-invasive medical devices. The wipes may be used in healthcare facilities, surgeries, nursing homes, and other areas, according to the intended use of the soaking agent.



Dry wipes



For small and large surfaces



Durable and convenient bucket



Non-dusting and lint-free white wipes



Recommended products

- ▶ Velox® Oxy ETA → p.6
- ▶ Velox® Top AF → p.7
- ▶ Velox® Spray → p.8
- ▶ Velox® Foam Extra → p.13
- ▶ Quatrodes® Forte → p.16
- ▶ Quatrodes® Extra → p.17
- ▶ Quatrodes® One → p.18
- ▶ Velodes® Silk → p.27

Packaging units:

Box with wipes for soaking insert-100 pcs
Insert for the box-100 pcs

Wipe grammage: 45 g/m²

Wipe size: 16 x 30 cm

Bucket material: HDPE with high-resistance to chemicals

Medical device class I



Additional information



Dosage:
0,5% → 5 ml
in 995 ml of water



Logistic details

Cap.	Product code	Scent	Pcs / Box	No of boxes / pallet	Pcs / europallet
1L	SSE-GOT-ML306	-	12	8	384
5L	SSE-GOT-ML039	-	1	30	120

Quatrodes® Forte

Concentrate for surfaces cleaning and disinfection

Product features:

- ▶ Concentration from 0.5% in 15 min
- ▶ Active against enveloped viruses (including Vaccinia, BVDV) Adeno, Polio. Bacteria, Fungi, Tuberculosis
- ▶ Recommended for ABS and acrylic dental and medical equipment, units, upholstery materials, countertops, floors, treatment beds
- ▶ High material tolerance on ABS plastic and upholstery materials - Famed opinion
- ▶ Does not require rinsing with water
- ▶ For disinfecting surfaces having contact with food
- ▶ Free from aldehyde and phenol
- ▶ 1 L of concentrate = 200 L of working solution
- ▶ Possible registrations: medical device and biocide

Directions for use:

Prepare a working solution by diluting the concentrate with water in a desired proportion. The choice of concentration depends on the desired spectrum and duration of operation. In order to obtain a concentration of 1% dilute 10 ml of the concentrate in 990 ml of water. Disinfected surfaces do not require rinsing with water. Do not combine with other detergents nor disinfectants. For professional use only.

Application:

Quatrodes® Forte is a liquid concentrate for simultaneous cleaning and disinfection of surfaces in the medical area (units, operating chairs and beds, operating tables, medical equipment, floors, walls), in gastronomy nursing homes, sanatoriums, gyms, fitness and SPA centres, hairdressing, massage and beauty salons and solariums (for disinfecting sunbeds), tattoos parlours, kindergartens and schools. The product is suitable for surfaces having contact with food. Product has a bactericidal, mycobactericidal, yeast- and virucidal effect. The substance can be used on all floor and above floor surfaces (doors, tabletops, cabinets, tables) and for cleaning and disinfecting hairdressing accessories.



Concentrate



Active against
B, F, V, Tbc,
Adeno, Polio



Medical and
non-medical
surfaces



Concentration
from 0.5%
in 15 min

Active substances in 100g of product:

3.76 g of N- (3-aminopropyl) -
N-dodecylpropane-1, 3-diamine, 3.39 g
of Poly (oxy-1,2-ethanedyl), alpha.
[2- (didecylmetyloamino) ethyl] -
omega.-hydroxy-, propanoate (salt),
0.75 g 2-aminoethanol

Effectiveness

	EN 14885	Time	
		15 min	30 min
Bacteria	EN 13727	0,5%*, 1%	
	EN 14561	1%	0,5%
	EN 13697	1%*	
Fungi (C. albicans)	EN 13624	0,5%	
	EN 14562	1%	0,5%
	EN 13697	1%*	
Mycobacterium tuberculosis (M. avium, M.terrae)	EN 14348	1%*,4%	1%
	EN 14563	4%	1%
Enveloped viruses	EN 14476	0,5%	
Adeno Virus	EN 14476	1%	
Polio Virus	EN 14476	4%	
BVDV Virus	EN 14476	0,5%	

* clean conditions



Matching products

Pumps

▶ MS-CHE-POMS → p.34

Packaging units: 1L bottle with doser, 5L canister

CE 2274 class IIa



Additional information



Dosage:
0,5% → 5 ml
in 995 ml of water



Matching products

Pumps

► MS-CHE-POMS → p.34

Logistic details

Cap.	Product code	Scent	Pcs / Box	No of boxes / pallet	Pcs europallet
1L	SSE-GOT-ML344	-	12	8	384
5L	SSE-GOT-ML351	-	1	30	120

Quatroles® Extra

Concentrate for surfaces cleaning and disinfection

Product features:

- Concentration from 0.5% in 15 min
- Active against Bacteria, Fungi, Tuberculosis, enveloped viruses (incl. Vaccinia, BVDV), Adeno, Polio
- Recommended for dental and medical equipment made of ABS and acrylic, units, beds and treatment chairs, tabletops, floors, operating tables
- High material tolerance to ABS and upholstery materials - Famed opinion
- Free from aldehydes and phenols, which doesn't discolour the disinfected surfaces
- For disinfecting surfaces having contact with food
- 1L of concentrate = 200 L of working solution
- Possible registrations: medical device and biocide

Directions for use:

Prepare a working solution by diluting the concentrate with water in a desired proportion. The choice of concentration depends on the desired spectrum and duration of operation. In order to obtain 0,5% concentration, dilute 5 ml of the concentrate in 995 ml of water. Disinfected surfaces do not require rinsing with water. When disinfected surfaces are to be made contact with food, rinse them with drinking water. Do not combine with other detergents nor disinfectants. For professional use only.

Application:

Quatroles® Extra is a liquid concentrate for simultaneous cleaning and disinfection of surfaces in the medical area (units, operating chairs and beds, operating tables, medical equipment, floors, walls), in gastronomy, nursing homes, sanatoriums, gyms, fitness and SPA centres, hairdressing, massage and beauty salons and solariums (for disinfecting sunbeds), tattoos parlours, kindergartens and schools. The product is suitable for surfaces having contact with food. Substance can be used on all floor and above floor surfaces (doors, tabletops, cabinets, tables) as well as for cleaning and disinfecting hairdressing accessories.



Concentrate



Active against B, F, V,
Tbc, Adeno, Polio,
Noro



Medical and
non-medical
surfaces



Concentration
from 0.5%
in 15 min

Active substances in 100g of product:

1.5 g N- (3-aminopropyl) -
N-dodecylpropane-1,3-diamine, 3.12 g
Poly (oxy-1,2-ethanediyl), alpha. [2-
(didecylmethylamino) ethyl] -omega-
hydroxy-, propanoate (salt), 1.25 g
Didecylmethylammonium chloride

Effectiveness

	EN 14885	Time		
		15 min	30 min	60 min
Bacteria	EN 13727	0,25%*, 0,5%		
	EN 14561	1%	0,5%	
	EN 16615	1%		
Bacteria (P. aeruginosa)	EN 13697	0,5%		
	EN 13624	0,25%*, 0,5%		
Fungi (C. albicans)	EN 14562	1%	0,5%	
	EN 13697	0,5%		
	EN 16615	1%		
Fungi (A. brasiliensis)	EN 13624			1,5%
Mycobacterium tuberculosis (M. terrae, M. avium)	EN 14348	0,5%*, 1%	0,5%	
	EN 14563	1%*, 2%	0,5%	
Enveloped viruses	EN 14476	0,5%		
Adeno Virus	EN 14476	0,5%*, 1%*	0,5%	
Polio Virus	EN 14476	1%	0,5%*	
Rota Virus	EN 14476	1%		
BVDV Virus	EN 14476	0,5%		
Noro Virus	EN 14476	4%		


* clean conditions



Packaging units: 1L bottle with doser, 5L canister

CE 2274 class IIa

Additional information

 Dosage:
1% → 10 ml
in 990 ml of water



Logistic details

Cap.	Product code	Scent	Pcs / Box	No of boxes / pallet	Pcs / europallet
1L	SSE-GOT-ML164	-	12	8	384
5L	SSE-GOT-ML171	-	1	30	120

Quatroles® Unit NF

Concentrate for cleaning and disinfection of suction devices

Product features:

- ▶ Concentration from 1% in 15 min
- ▶ Active against Bacteria, Fungi, Tuberculosis, enveloped viruses
- ▶ Recommended for suction systems, dental unit bowls, amalgam separators
- ▶ Does not foam
- ▶ For disinfection of a closed circuit in hydromassage tubs
- ▶ It prevents biofilm from forming

Directions for use:


The concentrate should be diluted with water. In order to obtain 2% concentration, dilute 20 ml of the product in 980 ml of water.


IN SUCTION SYSTEM: Introduce the solution to the suction system, depending on the degree of soiling, once or twice a day. Suck the whole product, turn the device and leave it for at least 15 minutes. After disinfection, flush the suction system with water. Rinsing can be performed after the recommended disinfection time or the next day, immediately after the unit is put into operation. Do not mix with other cleaning or disinfecting agents.


IN HYDROMASSAGE TUB: Before the first use, rinse the device with water to remove any possible residues of other product. Fill the tub with clean, cold water to the level above of the nozzles. Add 200 ml of the product for every 10 L of water. Turn the hydromassage system on for 15 minutes. After use, empty the tub and rinse with water. Do not combine with other agents. Do not use the product during the patient's healing bath.


Application:

The concentrate is designed for simultaneous cleaning and disinfection of suction systems, amalgam separators and dental unit bowls. The product does not foam. Used daily, it prevents biofilm from forming. It can also be used to disinfect closed circuit in hydromassage tubs that are medical devices. For professional use only.

Concentrate

Active against
B, F, V, Tbc

Dental unit

Concentration
from 0.5%
in 15 min

Active substances in 100g of product:

4.4 g dioctyl-dimethylammonium chloride, 0.6 g dodecyl-dimethylbenzylammonium chloride

Effectiveness	EN 14885	Time	
		15 min	60 min
Bacteria	EN 13727	0.5%*, 1%	
Fungi (C. albicans)	EN 13624	0,5%	
Mycobacterium tuberculosis (M. tuberculosis)	EN 14348		4%
Enveloped viruses	EN 14476	2%	

* clean conditions



Matching products

Pumps

▶ MS-CHE-POMS → p.34

Packaging units: 1L bottle with doser, 5L canister

CE 2274 class IIa

Additional information



Dosage:
0,5% → 5 ml
in 995 ml of water



Logistic details

Cap.	Product code	Scent	Pcs / Box	No of boxes / pallet	Pcs / europallet
1L	SSE-GOT-ML276	-	12	8	384
5L	SSE-GOT-ML283	-	1	30	120

Viruton® Forte

Concentrate for cleaning and disinfection of surgical instruments

Product features:

- ▶ Concentration from 0.5% in 15 min
- ▶ Active against Bacteria, Fungi, Tuberculosis, enveloped Viruses (incl. Vaccinia, BVDV), Adeno and Polio
- ▶ Anticorrosion properties and high material tolerance - Chirmed opinion (instruments manufacturer)
- ▶ For instruments made of steel, nickel, copper, aluminum
- ▶ Cleans and flushes blood, protein and secretions
- ▶ Recommended for endoscopes - Vimex opinion (endoscopes manufacturer)
- ▶ For ultrasonic cleaners

Directions for use:

Prepare working solution by diluting the concentrate with water. The level of dilution depends on a desired spectrum and duration of action. In order to obtain a concentration of 1% dilute 10 ml of the concentrate in 990 ml of water. Instruments shall be immersed in the solution immediately after use. All parts and holes must be completely filled in and covered by the solution. Once disinfected, rinse the instruments with water and dry. The uncontaminated working solution is active 14 days. For hygiene reasons, it is recommended to exchange the working solution on a daily basis. Do not combine with other detergents nor disinfectants. For professional use only. Read the product label before use.

Application:

An efficient concentrate for simultaneous cleaning and disinfection of surgical and rotary instruments. It has very good cleaning properties and a broad microbicidal efficacy in a short time of action. It is designed for manual disinfection and in ultrasonic washers. The product shows anticorrosive properties and high material tolerance. It is recommended for instruments made of stainless steel, nickel, copper, aluminum, rubber, porcelain, glass and plastic.



Concentrate



Active against B, F, V, Tbc, Adeno, Polio



Instruments



Concentration from 0.5% in 15 min

Active substances in 100g of product:

3.76 g of N-(3-aminopropyl) -N-dodecylpropane-1,3-diamine, 3.39 g of Poly (oxy- 1,2-ethanedyl), alpha. [2(didecylmethylamino) ethyl] - omega-hydroxy-, propanoate (salt), 0.75 g 2-aminoethanol

Effectiveness

	EN 14885	Time	
		15 min	30 min
Bacteria	EN 13727	0,5%*, 1%	
	EN 14561	1%	0,5%
	EN 13697	1%*	
Fungi (C. albicans)	EN 13624	0,5%	
	EN 14562	1%	0,5%
	EN 13697	1%*	
Mycobacterium tuberculosis (M.terrae, M. avium)	EN 14348	4%	1%
	EN 14563	4%	1%
Enveloped viruses	EN 14476	0,5%	
Adeno Virus	EN 14476	1%	
Polio Virus	EN 14476	4%	
BVDV Virus	EN 14476	0,5%	

* clean conditions



Packaging units: 1L bottle with doser, 5L canister

CE 2274 class IIB

Additional information



Dosage:
0,5% → 5 ml
in 995 ml of water



Logistic details

Cap.	Product code	Scent	Pcs / Box	No of boxes / pallet	Pcs / europallet
1L	SSE-43-GOT-ML382	-	12	8	384
5L	SSE-43-GOT-ML399	-	1	30	120

Viruton® Extra

Concentrate for cleaning and disinfection of surgical instruments

Product features:

- ▶ Concentration from 0.5% in 15 min
- ▶ Active against: Bacteria, Active against Fungi, Mycobacteria, enveloped viruses (incl. Vaccinia, BVDV), Adeno, Polio
- ▶ The protease enzyme ensures excellent cleaning properties
- ▶ For instruments made of steel, nickel, copper, aluminum
- ▶ For ultrasonic cleaners
- ▶ Recommended for endoscopes - Vimex opinion
- ▶ Anticorrosion properties and high material tolerance - Chirmed opinion
- ▶ 1 L concentrate = 200 L working solution

Directions for use:

Prepare a working solution by diluting the concentrate with water to a desired concentration. The level of dilution depends on a desired spectrum and duration of action. In order to obtain 0.5% concentration, dilute 5 ml of the concentrate with 995 ml of water. Immerse the instruments in the solution immediately after use. All parts and holes must be completely filled in covered by the solution. Once disinfected, reuse the instruments with water and dry. After disinfection rinse the instruments with water and dry them. The uncontaminated working solution is active for 14 days. For hygiene reasons, it is recommended to exchange the working solution on a daily basis. Do not combine with other detergents or disinfectants. For professional use only.

Application:

Viruton® Extra is a highly effective concentrate for simultaneous cleaning and disinfecting surgical and rotary instruments. Designed for manual disinfection, as well as in ultrasonic cleaners. It is recommended for instruments made of stainless steel, nickel, copper, aluminum, rubber, porcelain, glass and plastic.



Concentrate



Active against
Adeno, Polio



Contains protease
enzyme



Concentration
from 0.5%
in 15 min

Active substances in 100g of product:

11.5 g N- (3-aminopropyl) - N-
dodecylpropane-1,3-diamine, 3.12 g
Poly (oxy-1,2- ethanediyl), alpha. [2
(didecylmethylamino) ethyl] - omega-
hydroxy-, propanoate (salt), 1.25 g
Didecylmethylammonium chloride

Effectiveness

	EN 14885	Time		
		15 min	30 min	60 min
Bacteria	EN 13727	0,5%		
	EN 14561	1%	0,5%	
Fungi (C. albicans)	EN 13624	0,5%		
	EN 14562	1%	0,5%	
Fungi (A. brasiliensis)	EN 13624			1,5%
Mycobacterium tuberculosis (M.terrae, M. avium)	EN 14348	0,5%*, 1%	0,5%	
	EN 14563	1%*, 2%	0,5%	
Enveloped viruses	EN 14476	0,5%		
Adeno Virus	EN 14476	0,5%*	0,5%	
Polio Virus	EN 14476	1%	0,5%*	
BVDV Virus	EN 14476	0,5%		
Noro Virus	EN 14476	4%		2%

* clean conditions



Packaging units: 1L bottle with doser, 5L canister

CE 2274 class IIb



Dosage:
1% → 10 g
in 990 ml of water



Logistic details

Cap.	Product code	Scent	Pcs / Box	No of boxes / pallet	Pcs / europallet
1 kg	SSE-GOT-ML218	-	1	24	240
5 kg	SSE-GOT-ML225	-	1	16	80

Viruton® Pulver

Powder for cleaning and disinfecting tools, endoscopes, incubators and other surfaces of medical devices

Product features:

- ▶ Works already at a concentration of 1% in 10 min
- ▶ Active cleaning thanks to 4 enzymes
- ▶ Suitable for ultrasonic cleaners and endoscopes
- ▶ High-level disinfection (sporicidal activity)
- ▶ Powder formulation
- ▶ No activator required
- ▶ Exhibits anticorrosive properties and high material tolerance - opinion of the manufacturer of Chirmed tools

Directions for use:

According to the dosage table, measure the right amount of water at 25-30 °C, then pour in the right amount of powder using the provided measuring cup, stir the solution, cover the container and let it activate for 15 minutes. After this time, stir the solution again. The use of warm water (temp. 35 °C) will speed up the dissolution of the powder. Fine residues at the bottom of the container guarantee the stability of the solution over time. **Washing and disinfection of instruments and endoscopes:** thoroughly immerse the instruments observing the complete filling of all components and channels of medical equipment with the solution. After the specified time of action (10 min or 30 min depending on the chosen concentration), remove the tools from the solution and rinse with water of drinking water quality (preferably sterile distilled water or demineralized water). Change the solution daily. The product is not intended for copper and brass instruments, nor for mechanically damaged chrome and nickel-plated instruments. **Washing and disinfecting surfaces (depending on the chosen concentration):** spray the surfaces of medical equipment with the 1% solution or wipe with a cloth and leave to act for 10 min. The product leaves a white residue on the disinfected surface so it is recommended to wipe the residue with a disposable towel after disinfection. **Washing and disinfecting incubators:** detailed instructions for incubator disinfection available. Additional notes: Protect from moisture. Tightly close the container after use. Do not immerse wet instruments in powder. Do not use on disposable instruments. Follow the instructions of the manufacturer of the instruments and equipment. Do not use the product if the packaging is damaged. Store at temperatures below 25° C. Use test strips to check the activity of the solution over time.

Application:

Viruton® Pulver is a modern generation formulation for high-level cleaning and disinfection of invasive and non-invasive medical devices. It is recommended for surgical and dental instruments as well as all kinds of endoscopes, surgical probes, ventilator tubes and other anesthesia devices. It can also be used to disinfect the surfaces of medical equipment such as incubators. It is recommended for manual disinfection, in ultrasonic, automatic and in endoscope washers. Product for professional use.



Non - dusting



Active against Polio, Adeno, Noro, Spores



Active washing



Concentration from 1% in 10 min.

Active substances in 100g of product:

44 g sodium percarbonate, 26 g TAED, enzymes, organic acids, non-ionic surfactants, corrosion inhibitors.

Effectiveness

EN 14885

Time

10 min	15 min	30 min	60 min
	0,5%		
1%		0,5%	
		0,5%	
		0,5%	
1%		0,5%	
1%			0,5%
2%	1,5%	1%	
	0,5%		
			0,5%
1%			
	0,5%		
2%		1%	
2%		1%	
0,5%			
1%*, 2%		0,5%	
2%			
1%*		0,5%	
1%*, 2%			

* clean conditions



Packaging units: 1kg with measuring cup, 5kg with measuring cup (use the measuring cup according to the dosing table)

CE 2274 class IIB



Logistic details

Cap.	Product code	Scent	Pcs / Box	No of boxes / pallet	Pcs / europallet
1L	SSE-GOT-ML290	-	12	8	384
5L	SSE-GOT-ML337	-	1	30	120

Viruton® Bohr

Ready-to-use formula for cleaning and disinfection of rotary instruments

Product features:





- ▶ Active in 15 min. against Bacteria, Fungi, Tuberculosis, enveloped viruses (incl. Vaccinia, BVDV) Adeno and Polio
- ▶ Recommended for drills, bone cutters, abrasive and polishing tools, silicone rubbers, ceramic tiles and root canal treatment instruments
- ▶ Can be used by immersing in a tub or ultrasonic cleaner
- ▶ Contains inhibitors of corrosion
- ▶ For aluminum, steel, metal alloys
- ▶ Very good cleaning properties at dried organic residues

Directions for use:

Immerse the tools in undiluted preparation immediately after use. Carefully immerse all components. After the specified operating time, remove the tools from the bath or washcloth, rinse in water and dry immediately with a disposable napkin. In the event of visible contamination, remove them with a brush and immerse the disinfected product in the preparation again. After drying, sterilize clean instruments. Do not mix with other cleaning agents or disinfectants. The product is intended for use by a professional user

Application:

Formulation designed for the simultaneous cleaning and disinfection of precision rotary instruments and other invasive medical instruments. Viruton® Bohr is recommended for disinfection of drills, bone cutters made of steel, hard metals, diamonds, abrasive and polishing instruments such as silicone erasers, ceramic plates and instruments for root canal treatment. It can be used by immersion in a disinfection bath or in an ultrasonic cleaner.

	Liquid		Active against B, F, V, Tbc, Adeno, Polio		Rotary instruments		Active in 15 min.
Active substances in 100g of product:							
0.15 g N- (3-aminopropyl) -N-dodecylpropane-1,3-diamin, 0.14 g Poli (oxy-1,2-ethanedyl), alpha42- (didecylmethylamino) ethyl] - .omegahydroxy-, propanoate (salt)							
Effectiveness		EN 14885		High organic load			
Bacteria		EN 13727		15 min			
		EN 14561		15 min			
Fungi (C. albicans)		EN 13624		15 min			
		EN 14562		15 min			
Mycobacterium tuberculosis (M.terrae, M. avium)		EN 14348		15 min			
		EN 14563		15 min			
Enveloped viruses		EN 14476		15 min			
Adeno Virus		EN 14476		15 min			
Polio Virus		EN 14476		15 min			
BVDV Virus		EN 14476		15 min			



Packaging units: 1L bottle with doser, 5L canister

CE 2274 class IIb



Logistic details

Cap.	Product code	Scent	Pcs / Box	No of boxes / pallet	Pcs / europallet
1L	SSE-GOT-ML443	-	12	8	384

Viruton® Pre

Pre-disinfecting and cleaning foam for surgical instruments

Product features:

- ▶ Active in 15 min. against Bacteria, Fungi, Tuberculosis, enveloped viruses (incl. Vaccinia, BVDV), Adeno, Polio
- ▶ Designed for the pre- disinfection of instruments before the actual disinfection
- ▶ Recommended for medical and surgical instruments and containers
- ▶ Contains corrosion inhibitors
- ▶ High material tolerance
- ▶ Very good cleaning properties at dried organic residues thanks to a formula based on the synergy of surfactants and nanomicelle technology
- ▶ Preparation in the form of foam

Directions for use:

Place the tools immediately after use in the container and spray thoroughly with the undiluted product, making sure that all elements are covered with foam. After a certain period of time, remove the tools from the tray and rinse with water. Pre-disinfection process may take up to 48 hours. In case of visible contamination, clean the instruments with a brush and cover again with the foam. Once rinsed with water instruments shall be disinfected either manually or automatically. Do not combine with other detergents nor disinfectants. Use protective gloves.

Application:

The product is ideal for the initial disinfection and cleaning of instruments made of stainless steel, galvanized steel and aluminum, as well as materials made of plastics and rubber. It cleans well dried organic residues. Prevents dried organic soiling during collection and transport of instruments to the place of proper disinfection.



Foam



Active against
B, F, V, Tbc
Adeno, Polio



Tools



Active
in 15 min.

Active substances in 100g of product:

0,15 g N-(3- aminopropyl)-N-dodecylpropane 1,3-diamine, 0,14 g Poly(oxy-1,2- etanodilo),.alfa.-[2-(didecylmethyloamine) ethylo]-.omega.-hydroxy-,propanoate (salt)

Effectiveness

	EN 14885	Conditions	
		clean	dirty
Bacteria	EN 13727		15 min
	EN 14561		15 min
	EN 13697	15 min	
Fungi (C. albicans)	EN 13624		15 min
	EN 14562		15 min
	EN 13697	15 min	
Mycobacterium tuberculosis (M.terrae, M. avium)	EN 14348		15 min
	EN 14563		15 min
Enveloped viruses	EN 14476		15 min
Adeno Virus	EN 14476		15 min
Polio Virus	EN 14476		15 min
BVDV Virus	EN 14476		15 min



Packaging units: 1L spray bottle DUO

CE 2274 class IIB



Broad microcidal spectrum



NEW

Logistic details

Cap.	Product code	Scent	Pcs / Box	No of boxes / pallet	Pcs / europallet
500 ml	SSE-ALKGOT-ML516	-	24	8	1152
1L	SSE-ALKGOT-ML042	-	12	10	600
5L	SSE-ALKGOT-ML735	-	1	30	120

Velodes® Skin

Liquid for hygienic and surgical hand disinfection

Product features:

- Hygienic Hand Disinfection: 3 ml in 30 seconds
- Surgical Hand Disinfection: 2 x 3 ml in 2 x 60 seconds
- Extended surgical disinfection effect for up to 3 hours
- Thanks to additional protective substances, it protects the skin from drying
- Efficacy confirmed by tests according to the latest harmonized standards
- Reduction of the risk of infection in situations requiring a high sanitary regime

Directions for use:

HYGIENIC HAND DISINFECTION (according to EN 1500): rub 3 ml of the product in your hands for 30 seconds.

SURGICAL HAND DISINFECTION (according to EN 12791): Apply 3 ml of the product to clean, dry hands and rub in, apply further portions of 3 ml, so that the hands remain moist for a rubbing time of 120 seconds. The product also has a prolonged surgical disinfection effect up to 3 hours.

Application:

IN MEDICAL AREA

Velodes® Skin is recommended for hygienic and surgical hand disinfection in health care facilities, i.e. hospitals, clinics, treatment rooms, dental and others. Thanks to full virucidal activity, it is recommended for use in surgery, operating theaters and during surgical procedures in dentistry.

IN NON-MEDICAL AREA

In addition to the medical area, the product is recommended in all institutions and areas of life where it is advisable to maintain the highest standards of hygiene, such as hairdressing and beauty salons, nursing homes and fitness clubs, gyms, schools, offices, gastronomy, veterinary medicine.



Liquid



Hygienic disinfection in 90 sec.



Surgical disinfection in 90 sec.



Prolonged effect up to 3 hrs

Active substances in 100g of product:

80 g ethanol, 8 g propan-2-ol

Effectiveness

Bacteria

EN 14885

High organic load

Fungi (C. albicans)

EN 13727

15 sec.

EN 13624

30 sec.

EN 13697

3 min

Fungi (A. brasiliensis)

EN 13624

30 sec.

EN 13697

60 sec.

Mycobacterium tuberculosis (M. avium, M. terrae)

EN 14348

30 sec.

Enveloped viruses

EN 14476

15 sec.

Adeno Virus

EN 14476

15 sec.

Noro Virus

EN 14476

15 sec.

Polio Virus

EN 14476

30 sec.

Hygienic hand disinfection

EN 1500

3 ml in time 30 sec.

Surgical hand disinfection

EN 12791

2 x 3 ml in time 120 sec.
+ Prolonged effectiveness of surgical disinfection up to 3 hours



Matching products

Dispensers

- IM-DEZ-ML010 → p.33
- IM-DEZ-ML011 → p.33
- IM-DEZ-ML012 → p.33
- IM-DEZ-ML028P → p.33

- IM-DEZ-ML021D → p.34
- IM-DEZ-ML020D → p.34
- IM-DEZ-ML025D → p.34
- IM-DEZ-ML033D → p.34

- IM-DEZ-ML016 → p.34
- IN-CHG-0200 → p.34
- IN-CHG-0201 → p.34

Pumps

- IM-DEZ-ML557 → p.34
- IM-DEZ-ML558 → p.34
- IM-DEZ-ML559 → p.34
- IM-DEZ-ML560 → p.34

Packaging units: 500 ml bottle, 1L bottle, 1L refill, 5L canister



Logistic details

Cap.	Product code	Scent	Pcs / Box	No of boxes / pallet	Pcs / europallet
250 ml	SSE-43-GOT-ML535	-	40	8	1600
500 ml	SSE-43-GOT-ML542-24 MM	-	24	8	1152
1L	SSE-43-GOT-ML559	-	12	10	600
1L Refill	SSE-43-GOT-ML010	-	6	20	480
5L	SSE-43-GOT-ML566	-	1	30	120

Velodes Soft

Liquid for hygienic and surgical hand disinfection

Product features:

- ▶ Hygienic hand disinfection: 2 x 1,5 ml in 2x 15 sec.
- ▶ Surgical hand disinfection: 2 x 1,5 ml in 90 sec.
- ▶ Prolonged bactericidal effect up to 3 hours
- ▶ Active against Bacteria, Fungi, Mycobacteria, Viruses
- ▶ Skin-friendly pH
- ▶ Liquid formulation
- ▶ Dermatologically tested

Directions for use:


HYGIENIC HAND DISINFECTION (according to EN 1500): rub 2 x 1,5 ml of the product in your hands for 2 x 15 seconds.

SURGICAL HAND DISINFECTION (according to EN 12791): rub 2 x 1,5 ml of the product in clean and dry hands for 90 seconds. Keep the hands equally wet by adding more liquid. The product also has a prolonged surgical disinfection effect up to 3 hours.


For professional use only.

Application:


Velodes® Soft is a preparation, intended for hygienic and surgical hand disinfection with bactericidal, yeasticidal and limited virucidal activity. It works quickly and effectively against enveloped and non-enveloped viruses such as HIV, HBV, HCV, Rota virus.




Liquid



Hygienic hand disinfection in 30 sec.



Surgical hand disinfection in 90 sec.



Prolonged effect up to 3 hours

Active substances in 100g of product:

60 g propan-2-ol, 0.5 g chlorhexidine gluconate

Effectiveness	EN 14885	Conditions	
		clean	dirty
Bacteria	EN 13727	30 sec.	60 sec.
Fungi (C. albicans)	EN 1275	30 sec.	
	EN 13624	30 sec.	
Mycobacterium tuberculosis (M.terrae)	EN 14348	30 sec.	
Enveloped viruses	EN 14476	30 sec.	
Rota Virus	EN 14476	45 sec.	
BVDV Virus	EN 14476	30 sec.	
Hygienic hand disinfection	EN 1500	2 x 1,5 ml in 2 x 15 sec.	
Surgical hand disinfection	EN 12791	2 x 1,5 ml in 90 sec. + prolonged effect up to 3 hrs.	



Matching products

Dispensers

- ▶ IM-DEZ-ML009 → p.33
- ▶ IM-DEZ-ML007 → p.33
- ▶ IM-DEZ-ML010 → p.33
- ▶ IM-DEZ-ML011 → p.33
- ▶ IM-DEZ-ML012 → p.33
- ▶ IM-DEZ-ML028P → p.33
- ▶ IM-DEZ-ML021D → p.34

- ▶ IM-DEZ-ML020D → p.34
- ▶ IM-DEZ-ML025D → p.34
- ▶ IM-DEZ-ML033D → p.34
- ▶ IM-DEZ-ML016 → p.34
- ▶ IN-CHG-0200 → p.34
- ▶ IN-CHG-0201 → p.34

Pumps

- ▶ IM-DEZ-ML557 → p.34
- ▶ IM-DEZ-ML558 → p.34
- ▶ IM-DEZ-ML559 → p.34
- ▶ IM-DEZ-ML560 → p.34



Packaging units: 250 ml spray bottle, 500 ml bottle, 1L bottle, 1L refill, 5L canister



Logistic details

Cap.	Product code	Scent	Pcs / Box	No of boxes / pallet	Pcs / europallet
500 ml	SSE-ALKGOT-ML754-24 MM	-	24	8	1152
1L Refill	SSE-ALKGOT-ML942	-	6	20	480
5L	SSE-ALKGOT-ML080	-	1	30	120

Velodes Gel

Gel for hygienic and surgical hand disinfection

Product features:

- Hygienic hand disinfection: 2 x 1,5 ml in 30 sec. (2 x 15 sec.)
- Surgical hand disinfection: 2 x 1,5 ml in 2x 90 sec.
- Active against Rota and Noro viruses
- Skin-friendly pH
- Dermatologically tested gel
- Prolonged bactericidal effect up to 3 hours

Directions for use:

HYGIENIC HAND DISINFECTION (according to EN 1500): rub 2 x 1,5 ml of the product in your hands for 30 seconds. (2x15 sec.)

SURGICAL HAND DISINFECTION (according to EN 12791): rub 2 x 1,5 ml of the product in clean and dry hands for 90 seconds. Keep the hands equally wet by adding more liquid. The product also has a prolonged surgical disinfection effect up to 3 hours.

Application:

Velodes® Gel is a gel designed for hygienic and surgical hand disinfection. It has a bactericidal, yeastcidal, tuberculocidal and limited virucidal effect.

Active substances in 100g of product:	Effectiveness	EN 14885	Conditions	
			clean	dirty
63.7 g ethanol, 6.3 g propan-2-ol	Bacteria	EN 13727	30 sec.	30 sec.
	Fungi (C. albicans)	EN 13624	15 sec.	15 sec.
	Mycobacterium tuberculosis (M.terrae)	EN 14348	30 sec.	30 sec.
	Enveloped viruses	EN 14476	30 sec.	60 sec.
	BVDV Virus	EN 14476	30 sec.	30 sec.
	Rota Virus	EN 14476	30 sec.	60 sec.
	Noro Virus	EN 14476	30 sec.	30 sec.
	Hygienic hand disinfection	EN 1500	2 x 1,5 ml in 2 x 15 sec.	
	Surgical hand disinfection	EN 12791	2 x 1,5 ml in 2 x 90 sec. + prolonged effect up to 3 hrs.	



Matching products

Dispensers

- IM-DEZ-ML009 → p.33
- IM-DEZ-ML007 → p.33
- IM-DEZ-ML011 → p.33
- IM-DEZ-ML012 → p.33
- IM-DEZ-ML020D → p.34
- IM-DEZ-ML016 → p.34
- IN-CHG-0200 → p.34
- IN-CHG-0201 → p.34

Pumps

- IM-DEZ-ML558 → p.34
- IM-DEZ-ML560 → p.34



Packaging units: 500ml bottle, 1L bag, 5L canister

Available scents:



Händedesinfektion / Hand disinfection				
Hygienisch / Hygienic		Chirurgisch / Surgical		
Einwirkungszeit / Exposure time (min)		Einwirkungszeit / Exposure time (min)		
1/2	1	1-2,5	3	5
konz. 2 x 3 ml*				

*Keep hands well moistened with the concentrated product during the entire application time.

Logistic details

Cap.	Product code	Scent	Pcs / Box	No of boxes / pallet	Pcs / europallet
250 ml	SSE-ALKGOT-ML533	Tea tonic	40	8	1600
500 ml	SSE-ALKGOT-ML540		24	8	1152
1L	SSE-ALKGOT-ML557		12	10	600
1L Refill	SSE-ALKGOT-ML0105		6	20	480
5L	SSE-ALKGOT-ML564		1	30	120

Velodes Silk

Liquid for hygienic and surgical hand disinfection

Product features:

- ▶ Hygienic hand disinfection: 2 x 1,5 ml in 2 x 15 sec.
- ▶ Surgical hand disinfection: 2 x 1,5 ml in 2x 90 sec.
- ▶ Extended surgical disinfection effect for up to 3 hours
- ▶ Active against Rota, Noro
- ▶ Skin-friendly pH.
- ▶ Contains skin care substances
- ▶ Dermatologically tested
- ▶ Prevents the feeling of hands stickiness

Directions for use:

HYGIENIC HAND DISINFECTION (according to EN 1500): 2 x 1,5 ml of the product to be rubbed for 30 seconds (2 x 15 sec.).

SURGICAL HAND DISINFECTION (according to EN 12791): 2 x 1,5 ml of the product to be rubbed for 3 minutes (2 x 90 sec.). Keep hands equally wet by adding more liquid. The product also has a prolonged surgical disinfection effect up to 3 hours.

SURFACE DISINFECTION: Spray the surface from the distance of approximately 30 cm and leave it to dry. Do not apply on alcohol-sensitive surfaces (e.g. acrylic glass).

Application:

Ready to use, alcohol - based liquid for hygienic and surgical hand disinfection in the medical field, the food sector, industrial and domestic environments and also in public utilities. For surfaces having and not having contact with food. Designed also for disinfection of non-invasive medical devices.

Active substances in 100g of product:	Effectiveness	EN 14885	Conditions	
			clean	dirty
63.7 g ethanol, 6.3 g propan-2-ol	Bacteria	EN 13727	30 sec.	30 sec.
		EN 13697	30 sec.	60 sec.
	Fungi (C. albicans)	EN 13624	15 sec.	15 sec.
		EN 13697	30 sec.	
	Mycobacterium tuberculosis (M.terrae)	EN 14348	30 sec.	30 sec.
	Enveloped viruses	EN 14476	30 sec.	60 sec.
	BVDV Virus	EN 14476	30 sec.	30 sec.
	Rota Virus	EN 14476	30 sec.	60 sec.
	Noro Virus	EN 14476	30 sec.	30 sec.
	Hygienic hand disinfection	EN 1500	2 x 1,5 ml in 2 x 15 sec.	
	Surgical hand disinfection	EN 12791	2 x 1.5 ml in 2 x 90 sec. + prolonged effect up to 3 hrs.	
	Hygienic hand disinfection	VAH	2 x 3 ml in 2 x 30 sec.	



Matching products

Dispensers

- ▶ IM-DEZ-ML009 → p.33
- ▶ IM-DEZ-ML007 → p.33
- ▶ IM-DEZ-ML010 → p.33
- ▶ IM-DEZ-ML011 → p.33
- ▶ IM-DEZ-ML012 → p.33
- ▶ IM-DEZ-ML028P → p.33
- ▶ IM-DEZ-ML021D → p.34
- ▶ IM-DEZ-ML020D → p.34
- ▶ IM-DEZ-ML025D → p.34
- ▶ IM-DEZ-ML033D → p.34
- ▶ IM-DEZ-ML016 → p.34
- ▶ IN-CHG-0200 → p.34
- ▶ IN-CHG-0201 → p.34

Pumps

- ▶ IM-DEZ-ML557 → p.34
- ▶ IM-DEZ-ML558 → p.34
- ▶ IM-DEZ-ML559 → p.34
- ▶ IM-DEZ-ML560 → p.34



Packaging units: 250 ml spray bottle, 500 ml bottle, 1L bottle, 1L refill, 5L canister

CE 2274 class IIa



Logistic details

Cap.	Product code	Scent	Pcs / Box	No of boxes / pallet	Pcs / europallet
500 ml	SSE-GOT-ML092-500ML	-	24	8	1152
5L	SSE-GOT-ML078-5L	-	1	30	120

Velodes Scrub

Antibacterial emulsion for hygienic hand washing

Product features:

- ▶ Recommended for people with allergies for hand and body wash
- ▶ Perfect for doctors, dentists, nurses and people with sensitive skin, who wash their hands frequently
- ▶ Doesn't contain dyes or scented compositions
- ▶ Neutral pH for skin
- ▶ Contains properties unfavourable to bacterial growth - tested according to EN 1499 norm
- ▶ Contains panthenol which gives the skin sensation of softness
- ▶ Gives the skin a smooth feeling
- ▶ Recommended for food industry and healthcare institutions

Directions for use:

Apply a small amount of liquid on wet hands, rub to foam and rinse with water. Avoid contact with eyes.

Application:

Velodes Scrub is a gentle emulsion for washing hands and the whole body. It contains no dyes or fragrance. Ideal for doctors, dentists, nurses and people with sensitive skin who frequently wash their hands. The skin-neutral pH and the content of glycerin and coconut oil derivative allow daily hand washing without causing irritation. The emulsion can be used in the food industry, health care institutions, schools, kindergartens and restaurants. It has a skin-friendly pH and mild cleaning agents. Dermatologically tested. Tested in accordance with EN1499 Hygienic hand washing.



Emulsion



Hand cleanliness maintenance



Hand washing and care

Contains glycerine and coconut oil



Matching products

Dispensers

- ▶ IM-DEZ-ML011 → p.33
- ▶ IM-DEZ-ML012 → p.33
- ▶ IM-DEZ-ML020D → p.33
- ▶ IM-DEZ-ML016 → p.34
- ▶ IM-DOZ-0001 → p.34
- ▶ IM-DOZ-0001-POP → p.34

Pumps

- ▶ IM-DEZ-ML554 → p.34





Logistic details

Cap.	Product code	Scent	Pcs / Box	No of boxes / pallet	Pcs / europallet
500 ml	SSE-GOT-ML160-500ML	-	24	8	1152
5L	SSE-GOT-ML146-5L	-	1	30	120

Velodes Soap

Emulsion for hand and body washing

Product features:

- ▶ Intended for hand washing before surgical and after hygienic hand disinfection
- ▶ Suitable for washing head and the entire body
- ▶ Protects the skin from drying
- ▶ Recommended for medical personnel exposed to skin maceration related to frequent hand washing
- ▶ Skin-friendly pH
- ▶ Recommended for washing patients in the prevention of bedsores

Directions for use:

HAND WASHING: Apply a small amount of the emulsion on wet hands and spread until it foams, then rinse with water. Perform activities according to the WHO hand washing technique.

WASHING IN THE SHOWER: Spread the emulsion on wet skin and hair. Wash your body and rinse under running water.

WASHING IN THE BATH: Pour the emulsion into the bath filled with water. Wash your skin and hair, and then rinse under running water.

Application:

Velodes Soap is an emulsion for washing hands before surgical and after hygienic hand disinfection. It is suitable for washing the head and the entire body. Due to the content of conditioning ingredients such as glycerin and coconut oil, it is particularly recommended for medical personnel exposed by frequent washing to maceration of the skin. The emulsion is also recommended for washing patients in the prevention of bedsores. Dermatologically tested.



Liquid



Body cleanliness maintenance



Hand washing and care

Contains glycerin and coconut oil



Matching products

Dispensers

- ▶ IM-DEZ-ML011 → p.33
- ▶ IM-DEZ-ML012 → p.33
- ▶ IM-DEZ-ML020D → p.33
- ▶ IM-DEZ-ML016 → p.34
- ▶ IM-DOZ-0001 → p.34
- ▶ IM-DOZ-0001-POP → p.34

Pumps

- ▶ IM-DEZ-ML554 → p.34



Available scents:



Neutral
scent



White flowers
scent

Additional information



Cosmetic

Logistic details

Cap.	Product code	Scent	Pcs / Box	No of boxes / pallet	Pcs europallet
1L Refill	SSE-GOT-ML386-1L	neutral	6	20	480
1L Refill	SSE-GOT-ML379-1L	white flowers			

Velodes Soap Foam

Foam for hand-washing

Product features:

- ▶ Mild liquid applied in the form of foam for hygienic hand washing
- ▶ The liquid does not contain soap, so it can be used by people with allergies and intolerance to soap-based products
- ▶ The liquid creates a white, durable and thick foam gently covering the washed hands
- ▶ Due to the content of nourishing and moisturizing ingredients, it is recommended for medical personnel exposed to frequent maceration of the skin
- ▶ Dermatologically tested

Directions for use:

HAND WASHING: Apply a small amount of the foam on wet hands, then rinse with water.

WASHING IN THE SHOWER: Spread the foam on wet skin and hair. Wash your body and rinse under running water.

Application:

Velodes® Soap Foam is a gentle liquid applied in the form of foam for hygienic hand-washing. The product does not contain soap so it can be used by people with allergies and intolerance to soap-based products. The liquid creates white, durable and thick foam which gently covers washed hands. Thanks to skin-nourishing and moisturizing ingredients, it is particularly recommended for medical personnel susceptible to skin maceration related to frequent hand washing. Dermatologically tested.



Foam



Hand washing
and care



Cares and
moisturizes

Hygienic hand
washing



Neutral



White flowers



Matching products

Dispensers

- ▶ IM-DEZ-ML009 → p.33
- ▶ IM-DEZ-ML007 → p.33

Packaging units: 1L refill for the dispenser

pH 6±0,5



Logistic details

Cap.	Product code	Scent	Pcs / Box	No of boxes / pallet	Pcs / europallet
100 ml	MS-CHE-ML481-100ML	-	50	12	4800
500 ml	MS-CHE-ML548-500ML	-	30	6	1080

Velodes Cream

Protective hand cream

Product features:

- ▶ Gentle cream for hands and body skin care
- ▶ Protects, moisturizes and regenerates
- ▶ Nourishes dry, irritated and rough skin of the hands
- ▶ Recommended for everyday use after frequent hand washing and after hygienic disinfection
- ▶ Recommended for medical articular personnel and people exposed to drying and maceration of the skin as a result of frequent washing and wearing gloves
- ▶ Provides hands with silky softness
- ▶ Absorbs quickly and brings an immediate feeling of comfort
- ▶ Contains vitamin E, glycerine and olive oil

Directions for use:

Apply a small amount of cream on clean, dry hands. Gently rub into skin of your hands and nails. Repeat the process after each hand washing.

Application:

Velodes® Cream is a gentle cream for hands and body skin care. Protects, moisturizes and regenerates. It nourishes dry, irritated and rough skin. It is recommended for everyday use after frequent hand washing and after hygienic disinfection. It is particularly recommended for medical personnel and people exposed to drying and maceration of the skin as a result of frequent washing and wearing gloves. Provides silky hands softness. Quick absorption brings an immediate feeling of comfort. It does not contain parabens. Dermatologically tested.



Cream



Moisturizes, protects and nourishes



Contains olive oil

Contains vitamin E, glycerine



Matching products

Dispensers

▶ IM-DEZ-ML020D → p.34

Pumps

▶ IM-DEZ-ML554 → p.34

Packaging units: 100 ml tube, 500 ml bottle



Logistic details

Cap.	Product code	Scent	Pcs / Box	No of boxes / pallet	Pcs / europallet
100 ml	MS-CHE-ML027-100ML	-	50	12	4800
500 ml	MS-CHE-ML034-500ML	-	30	6	1080

Velodes Cream Aloe Vera

Moisturising hand cream

Product features:

- ▶ Perfectly moisturizes, nourishes and protects hand skin
- ▶ The rich formula, containing glycerine, aloe extract, pro-vitamin B5 and vitamin E, makes hand skin soft and silky
- ▶ It is recommended for everyday use after frequent hand washing and after hygienic disinfection
- ▶ Moisturizes, relaxes and reduces the feeling of roughness
- ▶ Particular recommended for doctors, medical personnel and people exposed to drying and maceration of the skin as a result of frequent washing and wearing gloves

Directions for use:

Apply a small amount of cream on clean, gently rub into skin of the hands and nails. Repeat the process after each hand washing.

Application:

Velodes® Cream Aloe Vera moisturizes, nourishes the skin of hands and protects it from the adverse influence of external factors. Rich formula, containing glycerine, aloe extract, pro-vitamin B5 and vitamin E, makes hand skin soft and silky.

It is recommended for everyday use after frequent hand washing and after hygienic disinfection. Especially recommended for medical personnel and people exposed to drying and maceration of the skin as a result of frequent washing and wearing gloves. Do not apply under diaper area for children 3 years of age or younger. Dermatologically tested.



Cream

Moisturizes
protects
and nourishesContains
aloe extractContains
vitamin E,
glycerine

Matching products

Dispensers

▶ IM-DEZ-ML020D → p.34

Pumps

▶ IM-DEZ-ML554 → p.34

Packaging units: 100 ml tube, 500 ml bottle

Elbow dispenser
for closed-system 1L
white



Product code:
IM-DEZ-ML009

Elbow dispenser
for closed-system 1L
white - violet



Product code:
IM-DEZ-ML009F

Dispenser
for closed system 1L
white



Product code:
IM-DEZ-ML007

Dispenser
for closed system 1L
white - violet



Product code:
IM-DEZ-ML007F

Dispenser
for closed system 1L
black



Product code:
IM-DEZ-ML007C

ABS elbow dispenser 1L



Product code:
IM-DEZ-ML010

ABS elbow dispenser 500 ml



Product code:
IM-DEZ-ML011

ABS elbow dispenser 500 ml
with tray



Product code:
IM-DEZ-ML012

Automatic, stand-alone E28 PLUS dispenser

Product features:

- ▶ Touch-free and automatic dispensing of liquid with possible adjustment of the dispensed dose
- ▶ RGB LED status illumination
- ▶ Eco-mode reducing electricity consumption
- ▶ The dispenser is designed and manufactured in Poland
- ▶ Made of powder coated steel
- ▶ 2-in-1 mains/battery power



Recommended products

- ▶ Velodes® Soft → p.25
- ▶ Velodes® Silk → p.27








Additional information







Availability
confirmation
needed

Technical data

Product code	Pojemność	Width	Height	Depth
IM-DEZ-ML028P	5L	280 mm	1205 mm	265 mm

Accessories for the PROCLEAN ECO dispenser				
<p>Stainless steel dispenser PROCLEAN ECO 1L</p>  <p>22,5 cm 6,5 cm 9,5 cm</p> <p>Product code: IM-DEZ-ML021D</p>	<p>Stainless steel dispenser PROCLEAN ECO 500 ml</p>  <p>19,8 cm 5,5 cm 7,7 cm</p> <p>Product code: IM-DEZ-ML020D</p>	<p>Adapter for the dispenser</p>  <p>23 cm 10,8 cm</p> <p>Product code: IM-DEZ-ML030D</p>	<p>Stand with H base</p>  <p>111,5 cm 31 cm 36 cm</p> <p>Product code: IM-DEZ-ML028D</p>	<p>Stand with O base</p>  <p>111,5 cm 35,5 cm</p> <p>Product code: IM-DEZ-ML031D</p>


Accessories for the PROCLEAN 7 wall dispenser:			
<p>PROCLEAN 2 dispenser white automatic wall mounted 1L</p>  <p>54,7 cm 24 cm 26,1 cm</p> <p>Product code: IM-DEZ-ML025D</p>	<p>PROCLEAN 7 dispenser white automatic wall mounted 1L</p>  <p>32,3 cm 8,5 cm 15,5 cm</p> <p>Product code: IM-DEZ-ML033D</p>	<p>Stand</p>  <p>130,8 cm 35,5 cm</p> <p>Product code: IM-DEZ-ML034D</p>	<p>Drip tray</p>  <p>0,16 cm 19,5 cm 15,5 cm</p> <p>Product code: IM-DEZ-ML035D</p>

<p>ABS and aluminum plastic dispenser 500 ml</p>  <p>27 cm 1,6 cm 8,1 cm</p> <p>Product code: IM-DEZ-ML016</p>	<p>Short pump gel / cream</p>  <p>24 mm</p> <p>Product code: IM-DEZ-ML554</p>	<p>Three nozzle spray pump for a 500 ml bottle</p>  <p>24 mm</p> <p>Product code: IM-DEZ-ML557</p>	<p>Gel or liquid pump for 500 ml bottle</p>  <p>24 mm</p> <p>Product code: IM-DEZ-ML558</p>
--	---	---	---

<p>Three-nozzle spray pump for a 1L bottle</p>  <p>28 mm</p> <p>Product code: IM-DEZ-ML559</p>	<p>Pump for gel or liquid into a 1L bottle</p>  <p>28 mm</p> <p>Product code: IM-DEZ-ML560</p>	<p>Canister pump</p>  <p>Product code: MS-CHE-POMS</p>	<p>Canister sprayer</p>  <p>Product code: IM-AKC-0033</p>
--	--	---	---

How to use a doser?

1. *Open the cap over the dispenser*



2. *Press the bottle.*
The liquid flows through the tube to the appropriate level of the dispenser



3. *Tilt the bottle as shown in the picture and pour the discharged amount of concentrate into the water*



How to prepare a working solution?

In order to obtain selected concentration, measure the appropriate amount of water in a tub or a bucket, then add the amount of disinfectant concentrate given in the table.

Concentration [%]	4%	2%	1%	0,5%	0,25%
Concentrate amount [ml]	40	20	10	5	2,5
Water amount [ml]	960	980	990	995	997,5

Concentrate efficiency (liters of solution)

Concentration [%]	Efficiency from 1l of concentrate
0,25%	400 liters
0,5%	200 liters
1%	100 liters
2%	50 liters
3%	33 liters
4%	25 liters

Partner:

FOR MORE INFORMATION PLEASE CONTACT OUR EXPORT DEPARTMENT

Export Department

tel.: +48 81 535 22 66

fax: +48 81 535 22 37

e-mail: export@medisept.pl

www.medisept.biz



More about our products at
www.medisept.biz



Manufacturer:
MEDISEPT Sp. z o.o.
ul. Ludwika Spiessa 4
20-270 Lublin, Poland

rev. 2022.03