

**MEDISEPT**

MEDISEPT Sp. z o.o.
Konopnica 159c
21-030 Motycz - lubelskie - Polska
Phone: +48 81 535 22 92
g.gromadzki@medisept.pl
<https://medisept.pl/>

Ingredient data sheet

According to 648/2004/EC , 907/2006/EC

Product identifier: Quatrosdes Unit NF

10 % or more			
CAS	Common chemical name or IUPAC name	INCI	European Pharmacopoeia
7732-18-5	Water	AQUA	

1 % or over, but less than 10 %			
CAS	Common chemical name or IUPAC name	INCI	European Pharmacopoeia
34590-94-8	Dipropylene Glycol Methyl Ether	PPG-2 METHYL ETHER	
7320-34-5	Tetrapotassium pyrophosphate	TETRAPOTASSIUM PYROPHOSPHATE	
5538-94-3	Dimethyldioctylammonium chloride		
64-17-5	ethanol	ALCOHOL	
164462-16-2	Alanine, N,N-bis(carboxymethyl)-, trisodiumsalt	TRISODIUM DICARBOXYMETHYL ALANINATE	
1310-58-3	potassium hydroxide	POTASSIUM HYDROXIDE	

0,1 % or over, but less than 1 %			
CAS	Common chemical name or IUPAC name	INCI	European Pharmacopoeia
85409-22-9	Quaternary ammonium compounds, benzyl-C12-14-alkyldimethyl, chlorides	BENZALKONIUM CHLORIDE	