

**MEDISEPT**

MEDISEPT Sp. z o.o.
Konopnica 159c
21-030 Motycz - lubelskie - Polska
Phone: +48 81 535 22 92
g.gromadzki@medisept.pl
<https://medisept.pl/>

Ingredient data sheet

According to 648/2004/EC , 907/2006/EC

Product identifier: Velox Oxy ETA

| 10 % or more | | | |
|--------------|------------------------------------|---------|------------------------|
| CAS | Common chemical name or IUPAC name | INCI | European Pharmacopoeia |
| 64-17-5 | ethanol | ALCOHOL | |
| 7732-18-5 | Water | AQUA | |

| 1 % or over, but less than 10 % | | | |
|---------------------------------|------------------------------------|-------------------|------------------------|
| CAS | Common chemical name or IUPAC name | INCI | European Pharmacopoeia |
| 67-63-0 | propan-2-ol | ISOPROPYL ALCOHOL | |
| 7722-84-1 | Hydrogen peroxide solution | HYDROGEN PEROXIDE | |

| 0,1 % or over, but less than 1 % | | | |
|----------------------------------|--|-----------------------------------|------------------------|
| CAS | Common chemical name or IUPAC name | INCI | European Pharmacopoeia |
| 6419-19-8 | Nitritotrimethylenetris(phosphonic acid) | AMINOTRIMETHYLENE PHOSPHONIC ACID | |
| 54549-24-5 | Hexyl D-glucoside | | |
| 7664-38-2 | Phosphoric acid | PHOSPHORIC ACID | |

| less than 0,1 % | | | |
|-----------------|------------------------------------|------|------------------------|
| CAS | Common chemical name or IUPAC name | INCI | European Pharmacopoeia |
| 13598-36-2 | phosphonic acid | | |