

**MEDISEPT**

MEDISEPT Sp. z o.o.
Konopnica 159c
21-030 Motycz - lubelskie - Polska
Phone: +48 81 535 22 92
g.gromadzki@medisept.pl
<https://medisept.pl/>

Ingredient data sheet

According to 648/2004/EC , 907/2006/EC

Product identifier: Velox Top AF Grapefruit

10 % or more			
CAS	Common chemical name or IUPAC name	INCI	European Pharmacopoeia
67-63-0	propan-2-ol	ISOPROPYL ALCOHOL	
7732-18-5	Water	AQUA	
64-17-5	ethanol	ALCOHOL	

0,1 % or over, but less than 1 %			
CAS	Common chemical name or IUPAC name	INCI	European Pharmacopoeia
139734-65-9	Amines, n-C10-16-alkyltrimethylenedi-, reaction products with chloroacetic acid		
7173-51-5	Didecyldimethylammonium chloride	DIDECYLDIMONIUM CHLORIDE	
57-55-6	Propane-1,2-diol	PROPYLENE GLYCOL	

less than 0,1 %			
CAS	Common chemical name or IUPAC name	INCI	European Pharmacopoeia
5989-27-5	d-limonene	LIMONENE	
78-70-6	Linalool	LINALOOL	
	perfumes		