

**MEDISEPT**

MEDISEPT Sp. z o.o.  
Konopnica 159c  
21-030 Motycz - lubelskie - Polska  
Phone: +48 81 535 22 92  
g.gromadzki@medisept.pl  
<https://medisept.pl/>

## Ingredient data sheet

According to 648/2004/EC , 907/2006/EC

Product identifier: Viruton Pulver

10 % or more			
CAS	Common chemical name or IUPAC name	INCI	European Pharmacopoeia
15630-89-4	Disodium carbonate, compound with hydrogen peroxide (2:3)	SODIUM CARBONATE PEROXIDE	
10543-57-4	N,N'-ethylenebis[N-acetylacetylamide]	TAED	
124-04-9	adipic acid	ADIPIC ACID	

1 % or over, but less than 10 %			
CAS	Common chemical name or IUPAC name	INCI	European Pharmacopoeia
68131-40-8	Alcohols, C11-15-secondary, ethoxylated	C11-15 SEC-PARETH-12	
3794-83-0	tetrasodium (1-hydroxyethylidene)bisphosphonate	TETRASODIUM ETIDRONATE	
497-19-8	sodium carbonate	SODIUM CARBONATE	
7601-54-9	Trisodium orthophosphate	TRISODIUM PHOSPHATE	

0,1 % or over, but less than 1 %			
CAS	Common chemical name or IUPAC name	INCI	European Pharmacopoeia
10213-79-3	Disodium metasilicate· 5 H <sub>2</sub> O	SODIUM METASILICATE	
7732-18-5	Water	AQUA	

less than 0,1 %			
CAS	Common chemical name or IUPAC name	INCI	European Pharmacopoeia
25322-68-3	Polyethylenglycol 200		
	Vinasses, residue of fermentation containing biomass of bakers yeast ( <i>Saccharomyces cerevisiae</i> )		
9014-01-1	Subtilisin	SUBTILISIN	
9001-62-1	Lipase, triacylglycerol	LIPASE	
9000-90-2	Amylase, α-		